

KV β (H-300): sc-134864

BACKGROUND

Voltage-gated K⁺ channels in the plasma membrane control the repolarization and the frequency of action potentials in neurons, muscles and other excitable cells. The KV gene family encodes more than 30 proteins that comprise the subunits of the K⁺ channels, and they vary in their gating and permeation properties, subcellular distribution and expression patterns. Functional KV channels assemble as tetramers consisting of pore-forming α subunits (KV), which include the KV1, KV2, KV3 and KV4 proteins, and accessory or KV-subunits that modify the gating properties of the coexpressed KV subunits. KV β , also known as KCNAB1 (potassium voltage-gated channel, shaker-related subfamily, β member 1), is a 419 amino acid accessory K⁺ channel protein that exists as three alternatively spliced isoforms and regulates the activity of the pore-forming α subunit. It is expressed in the brain, with highest levels detected in the caudate nucleus, hippocampus and thalamus.

REFERENCES

1. Majumder, K., De Biasi, M., Wang, Z. and Wible, B.A. 1995. Molecular cloning and functional expression of a novel potassium channel β subunit from human atrium. *FEBS Lett.* 361: 13-16.
2. Morales, M.J., Castellino, R.C., Crews, A.L., Rasmusson, R.L. and Strauss, H.C. 1995. A novel β subunit increase potassium channel α subunits. *J. Biol. Chem.* 270: 6272-6277.
3. England, S.K., Uebele, V.N., Shear, H., Kodali, J., Bennett, P.B. and Tamkun, M.M. 1995. Characterization of a voltage-gated K⁺ channel β subunit expressed in human heart. *Proc. Natl. Acad. Sci. USA* 92: 6309-6313.
4. McCormack, K., McCormack, T., Tanouye, M., Rudy, B. and Stühmer, W. 1995. Alternative splicing of the human expression of the β 2 gene product. *FEBS Lett.* 370: 32-36.
5. England, S.K., Uebele, V.N., Kodali, J., Bennett, P.B. and Tamkun, M.M. 1996. A novel K⁺ channel β -subunit (hKV β 1.3) is produced via alternative mRNA splicing. *J. Biol. Chem.* 270: 28531-28534.
6. Leicher, T., Roeper, J., Weber, K., Wang, X. and Pongs, O. 1997. Structural and functional characterization of human potassium channel subunit β 1 (KCNA1B). *Neuropharmacology* 35: 787-795.
7. Leicher, T., Bähring, R., Isbrandt, D. and Pongs, O. 1999. Coexpression of the KCNA3B gene product with KV1.5 leads to a novel A-type potassium channel. *J. Biol. Chem.* 273: 35095-35101.
8. Manganas, L.N. and Trimmer, J.S. 2000. Subunit composition determines KV1 potassium channel surface expression. *J. Biol. Chem.* 275: 29685-29693.

SOURCE

KV β (H-300) is a rabbit polyclonal antibody raised against amino acids 120-419 mapping at the C-terminus of KV β .1 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

KV β (H-300) is recommended for detection of KV β .1, KV β .2, and to a lesser extent, KV β .3 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

KV β (H-300) is also recommended for detection of KV β .1, KV β .2, and to a lesser extent, KV β .3 in additional species, including equine, canine, bovine, porcine and avian.

Molecular Weight of KV β : 47 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker[™] compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker[™] Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz[™] Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.