

# AMID (H-300): sc-135000

## BACKGROUND

AMID (apoptosis-inducing factor (AIF)-like mitochondrion-associated inducer of death), also called p53-responsive gene 3 (PRG3), is a member of the FAD-dependent oxidoreductase family. AMID is a caspase-independent proapoptotic flavoprotein with NAD(P)H oxidase activity localizing to the cytosol and associated with the outer mitochondrial membrane. AMID shares significant homology with AIF and NADH-oxidoreductases. It is expressed in most normal tissues and its expression is upregulated by p53. Two AMID isoforms exist due to alternative splicing. Isoform 1 is the full length protein and isoform 2 is missing amino acids 99-138. Isoform 2 also has an additional three amino acids inserted after residue 206. Overexpression of AMID leads to apoptosis.

## REFERENCES

1. Wu, M., Xu, L.G., Li, X., Zhai, Z. and Shu, H.B. 2002. AMID, an apoptosis-inducing factor-homologous mitochondrion-associated protein, induces caspase-independent apoptosis. *J. Biol. Chem.* 277: 25617-25623.
2. Lü, C.X., Fan, T.J., Hu, G.B. and Cong, R.S. 2003. Apoptosis-inducing factor and apoptosis. *Sheng Wu Hua Xue Yu Sheng Wu Wu Li Xue Bao.* 35: 881-885.
3. Jan, Y., Matter, M., Pai, J.T., Chen, Y.L., Pilch, J., Komatsu, M., Ong, E., Fukuda, M. and Ruoslahti, E. 2004. A mitochondrial protein, Bit1, mediates apoptosis regulated by integrins and Groucho/TLE corepressors. *Cell* 116: 751-762.
4. Wu, M., Xu, L.G., Su, T., Tian, Y., Zhai, Z. and Shu, H.B. 2004. AMID is a p53-inducible gene downregulated in tumors. *Oncogene* 23: 6815-6819.
5. Marshall, K.R., Gong, M., Wodke, L., Lamb, J.H., Jones, D.J., Farmer, P.B., Scrutton, N.S. and Munro, A.W. 2005. The human apoptosis-inducing protein AMID is an oxidoreductase with a modified flavin cofactor and DNA binding activity. *J. Biol. Chem.* 280: 30735-30740.
6. Hu, Q., Chang, J., Tao, L., Yan, G., Xie, M. and Wang, Z. 2005. Endoplasmic reticulum mediated necrosis-like apoptosis of HeLa cells induced by Ca<sup>2+</sup> oscillation. *J. Biochem. Mol. Biol.* 38: 709-716.
7. Li, W., Sun, L., Liang, Q., Wang, J., Mo, W. and Zhou, B. 2006. Yeast AMID homologue Ndi1p displays respiration-restricted apoptotic activity and is involved in chronological aging. *Mol. Biol. Cell* 17: 1802-1811.
8. Mei, J., Webb, S., Zhang, B. and Shu, H.B. 2006. The p53-inducible apoptotic protein AMID is not required for normal development and tumor suppression. *Oncogene* 25: 849-856.
9. Varecha, M., Amrichová, J., Zimmermann, M., Ulman, V., Lukášová, E. and Kozubek, M. 2007. Bioinformatic and image analyses of the cellular localization of the apoptotic proteins endonuclease G, AIF, and AMID during apoptosis in human cells. *Apoptosis* 12: 1155-1171.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## CHROMOSOMAL LOCATION

Genetic locus: AIFM2 (human) mapping to 10q22.1; Aifm2 (mouse) mapping to 10 B4.

## SOURCE

AMID (H-300) is a rabbit polyclonal antibody raised against amino acids 1-300 mapping at the N-terminus of AMID of human origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

## APPLICATIONS

AMID (H-300) is recommended for detection of AIF-like mitochondrion-associated inducer of death of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for AMID siRNA (h): sc-72339, AMID siRNA (m): sc-72340, AMID shRNA Plasmid (h): sc-72339-SH, AMID shRNA Plasmid (m): sc-72340-SH, AMID shRNA (h) Lentiviral Particles: sc-72339-V and AMID shRNA (m) Lentiviral Particles: sc-72340-V.

Molecular Weight of AMID: 41 kDa.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.