

Latrophilin-2 (H-45): sc-135155

BACKGROUND

The latrophilin family of G protein-coupled receptors consists of three members, Latrophilin-1, 2 and 3, each of which displays distinct tissue distribution and function. Latrophilin-1, the most characterized member of this family, acts as a receptor for α -latrotoxin, a component of venom from the black widow spider. Binding of α -latrotoxin to Latrophilin-1 triggers synaptic vesicle exocytosis via both Ca^{2+} -dependent and -independent mechanisms, resulting in vesicle depletion. Latrophilin-1 is abundant in brain and present in endocrine cells. Latrophilin-3 is also brain-specific, whereas Latrophilin-2 expression is ubiquitous.

REFERENCES

1. Matsushita, H., et al. 1999. The latrophilin family: multiply spliced G protein-coupled receptors with differential tissue distribution. *FEBS Lett.* 443: 348-352.
2. Bittner, M.A., et al. 2000. α -latrotoxin and its receptors C1RL (latrophilin) and neurexin 1 α mediate effects on secretion through multiple mechanisms. *Biochimie* 82: 447-452.
3. Van Renterghem, C., et al. 2000. α -latrotoxin forms calcium-permeable membrane pores via interactions with latrophilin or neurexin. *Eur. J. Neurosci.* 12: 3953-3962.
4. Sudhof, T.C., et al. 2001. α -latrotoxin and its receptors: neurexins and C1RL/Latrophilins. *Annu. Rev. Neurosci.* 24: 933-962.
5. Nicholson, G.M., et al. 2002. Spiders of medical importance in the Asia-Pacific: atracotoxin, latrotoxin and related spider neurotoxins. *Clin. Exp. Pharmacol. Physiol.* 29: 785-794.
6. Ushkaryov, Y.A., et al. 2004. The multiple actions of black widow spider toxins and their selective use in neurosecretion studies. *Toxicon* 43: 527-542.

CHROMOSOMAL LOCATION

Genetic locus: LPHN2 (human) mapping to 1p31.1; Lphn2 (mouse) mapping to 3 H3.

SOURCE

Latrophilin-2 (H-45) is a rabbit polyclonal antibody raised against amino acids 326-370 mapping within an internal region of Latrophilin-2 of human origin.

PRODUCT

Each vial contains 200 μg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

APPLICATIONS

Latrophilin-2 (H-45) is recommended for detection of Latrophilin-2 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μg per 100-500 μg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Latrophilin-2 (H-45) is also recommended for detection of Latrophilin-2 in additional species, including equine, canine, bovine, porcine and avian.

Suitable for use as control antibody for Latrophilin-2 siRNA (h): sc-60919, Latrophilin-2 siRNA (m): sc-60920, Latrophilin-2 shRNA Plasmid (h): sc-60919-SH, Latrophilin-2 shRNA Plasmid (m): sc-60920-SH, Latrophilin-2 shRNA (h) Lentiviral Particles: sc-60919-V and Latrophilin-2 shRNA (m) Lentiviral Particles: sc-60920-V.


Molecular Weight of Latrophilin-2: 163 kDa.

Positive Controls: human smooth muscle extract: sc-363778.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



Latrophilin-2 (H-45): sc-135155. Western blot analysis of Latrophilin-2 expression in human smooth muscle tissue extract.

RESEARCH USE

For research use only, not for use in diagnostic procedures.



Try **Latrophilin-2 (E-3): sc-514197**, our highly recommended monoclonal alternative to Latrophilin-2 (H-45).