

DIAPH3 (H-42): sc-135261

BACKGROUND

DIAPH3 (diaphanous homolog 3), also known as DIAP3 or DRF3, is a 1,193 amino acid member of the formin homology protein family and is required for the correct function of various cellular processes. DIAPH3 binds to both profilin, a protein involved in cell maintenance, and to the GTP-bound form of Rho (Rho-GTP). Binding to both of these proteins allows DIAPH3 to recruit profilin to the membrane, in a Rho-dependent manner. At the membrane, DIAPH3 promotes actin polymerization and is required for stress fiber formation, cytokinesis and transcriptional activation of the serum response factor (SRF). DIAPH3 also regulates actin dynamics by coupling Src tyrosine kinase (c-Src) and Rho during actin signaling events. DIAPH3 contains one diaphanous autoregulatory domain (DAD) and one Rho GTPase-binding domain (GBD). When DAD and GBD are intramolecularly bound, the GBD is occupied and DIAPH3 is inactive. Interruption of the DAD-GBD bond allows the GBD to bind to Rho-GTP, thus activating DIAPH3. Seven isoforms of DIAPH3 exist due to alternative splicing events.

REFERENCES

- Peng, J., Wallar, B.J., Flanders, A., Swiatek, P.J. and Alberts, A.S. 2003. Disruption of the diaphanous-related formin Drf1 gene encoding mDia1 reveals a role for Drf3 as an effector for Cdc42. *Curr. Biol.* 13: 534-545.
- Katoh, M. and Katoh, M. 2004. Identification and characterization of human DIAPH3 gene in silico. *Int. J. Mol. Med.* 13: 473-478.
- Katoh, M. and Katoh, M. 2004. Identification and characterization of human FHOD3 gene in silico. *Int. J. Mol. Med.* 13: 615-620.
- Katoh, M. and Katoh, M. 2004. Identification and characterization of human FHDC1, mouse Fhdc1 and zebrafish fhdc1 genes in silico. *Int. J. Mol. Med.* 13: 929-934.
- Katoh, Y. and Katoh, M. 2004. Identification and characterization of CDC50A, CDC50B and CDC50C genes in silico. *Oncol. Rep.* 12: 939-943.

CHROMOSOMAL LOCATION

Genetic locus: DIAPH3 (human) mapping to 13q21.2; Diaph3 (mouse) mapping to 14 E1.

SOURCE

DIAPH3 (H-42) is a rabbit polyclonal antibody raised against amino acids 54-95 mapping near the N-terminus of DIAPH3 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

APPLICATIONS

DIAPH3 (H-42) is recommended for detection of DIAPH3 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

DIAPH3 (H-42) is also recommended for detection of DIAPH3 in additional species, including canine and bovine.


Suitable for use as control antibody for DIAPH3 siRNA (h): sc-62210, DIAPH3 siRNA (m): sc-62211, DIAPH3 shRNA Plasmid (h): sc-62210-SH, DIAPH3 shRNA Plasmid (m): sc-62211-SH, DIAPH3 shRNA (h) Lentiviral Particles: sc-62210-V and DIAPH3 shRNA (m) Lentiviral Particles: sc-62211-V.

Molecular Weight of DIAPH3: 137 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



DIAPH3 (H-42): sc-135261. Western blot analysis of DIAPH3 expression in non-transfected CHO (A) and mouse DIAPH3 transfected CHO (B) whole cell lysates.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

MONOS
Satisfaction
Guaranteed

Try **DIAPH3 (4D5): sc-293288**, our highly recommended monoclonal alternative to DIAPH3 (H-42).