

cep-1 (cC-18): sc-135460

The Power to Question

BACKGROUND

p53, a DNA-binding, oligomerization domain- and transcription activation domain-containing tumor suppressor, upregulates growth arrest and apoptosis-related genes in response to stress signals, thereby influencing programmed cell death, cell differentiation and cell cycle control mechanisms. cep-1 (p53-like protein), is the 644 amino acid *C. elegans* homolog of mammalian tumor suppressor p53. While the DNA-binding domains of p53 and cep-1 share little homology, both recognize nearly identical DNA sequences. Ubiquitously expressed in embryos, cep-1 is required to activate apoptosis of germ cells in response to genotoxic stress, and additionally activates *C. elegans* pgh-1, the homolog of Gas1 (growth arrest-specific 1). Together, cep-1 and pgh-1 mediate embryonic viability of mutations of TELO2, which is an essential gene in the regulation of checkpoint responses. Overexpression of cep-1 has been found to cause widespread caspase-independent cell death.

REFERENCES

1. Derry, W.B., Putzke, A.P. and Rothman, J.H. 2001. *Caenorhabditis elegans* p53: role in apoptosis, meiosis, and stress resistance. *Science* 294: 591-595.
2. Huyen, Y., Jeffrey, P.D., Derry, W.B., Rothman, J.H., Pavletich, N.P., Stavridi, E.S. and Halazonetis, T.D. 2004. Structural differences in the DNA binding domains of human p53 and its *C. elegans* ortholog cep-1. *Structure* 12: 1237-1243.
3. Schumacher, B., Hanazawa, M., Lee, M.H., Nayak, S., Volkmann, K., Hofmann, E.R., Hofmann, R., Hengartner, M., Schedl, T. and Gartner, A. 2005. Translational repression of *C. elegans* p53 by GLD-1 regulates DNA damage-induced apoptosis. *Cell* 120: 357-368.
4. D'Erchia, A.M., Tullo, A., Lefkimiatis, K., Saccone, C. and Sbisà, E. 2006. The fatty acid synthase gene is a conserved p53 family target from worm to human. *Cell Cycle* 5: 750-758.
5. Derry, W.B., Bierings, R., van Iersel, M., Satkunendran, T., Reinke, V. and Rothman, J.H. 2007. Regulation of developmental rate and germ cell proliferation in *Caenorhabditis elegans* by the p53 gene network. *Cell Death Differ.* 14: 662-670.
6. Stergiou, L., Doukoumetzidis, K., Sendoel, A. and Hengartner, M.O. 2007. The nucleotide excision repair pathway is required for UV-C-induced apoptosis in *Caenorhabditis elegans*. *Cell Death Differ.* 14: 1129-1138.
7. Arum, O. and Johnson, T.E. 2007. Reduced expression of the *Caenorhabditis elegans* p53 ortholog cep-1 results in increased longevity. *J. Gerontol. A Biol. Sci. Med. Sci.* 62: 951-959.
8. Tavernarakis, N., Pasparaki, A., Tasdemir, E., Maiuri, M.C. and Kroemer, G. 2008. The effects of p53 on whole organism longevity are mediated by autophagy. *Autophagy* 4: 870-873.
9. Gao, M.X., Liao, E.H., Yu, B., Wang, Y., Zhen, M. and Derry, W.B. 2008. The SCF FSN-1 ubiquitin ligase controls germline apoptosis through cep-1/p53 in *C. elegans*. *Cell Death Differ.* 15: 1054-1062.

SOURCE

cep-1 (cC-18) is an affinity purified rabbit polyclonal antibody raised against a peptide mapping near the C-terminus of cep-1 of *C. elegans* origin.

PRODUCT

Each vial contains 100 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-135460 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

cep-1 (cC-18) is recommended for detection of cep-1 of *Caenorhabditis elegans* origin by Western Blotting (starting dilution 1:100, dilution range 1:50-1:500), immunofluorescence (starting dilution 1:25, dilution range 1:25-1:250) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Molecular Weight of cep-1: 75 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.