

C/EBP δ (92.69): sc-135733

BACKGROUND

The transcription factor C/EBP α (CCAAT-enhancer binding protein) is a heat-stable, sequence-specific DNA-binding protein first purified from rat liver nuclei that binds avidly to several different *cis*-regulatory DNA sequences commonly associated with viral and cellular genes transcribed by RNA polymerase II. C/EBP α regulates gene expression in a variety of tissues including liver, adipose, lung and intestine. C/EBP α uses a bipartite structural motif to bind DNA. Two protein chains dimerize through a set of amphipathic α helices termed the leucine zipper. Highly basic polypeptide regions emerge from the zipper to form a linked set of DNA contact surfaces. C/EBP α appears to function exclusively in terminally-differentiated, growth-arrested cells. Additional family members include C/EBP β , C/EBP γ , C/EBP δ and C/EBP ϵ , all of which exhibit similar DNA-binding specificities and affinities to C/EBP α . Furthermore, C/EBP β and C/EBP δ readily form heterodimers both with each other as well as with C/EBP α .

REFERENCES

- Johnson, P.F., et al. 1987. Identification of a rat liver nuclear protein that binds to the enhancer core element of three animal viruses. *Genes Dev.* 1: 133-146.
- Landschulz, W.H., et al. 1988. Isolation of a recombinant copy of the gene encoding C/EBP. *Genes Dev.* 2: 786-800.
- Birkenmeier, E.H., et al. 1989. Tissue-specific expression, developmental regulation, and genetic mapping of the gene encoding CCAAT/enhancer binding protein. *Genes Dev.* 3: 1146-1156.
- Cao, Z., et al. 1991. Regulated expression of three C/EBP isoforms during adipose conversion of 3T3-L1 cells. *Genes Dev.* 5: 1538-1552.
- Umek, R.M., et al. 1991. CCAAT-enhancer binding protein: a component of a differentiation switch. *Science* 251: 288-292.
- Williams, S.C., et al. 1991. A family of C/EBP-related proteins capable of forming covalently linked leucine zipper dimers *in vitro*. *Genes Dev.* 5: 1553-1567.
- Davydov, I.V., et al. 1995. Cloning of the cDNA encoding human C/EBP γ , a protein binding to the PRE-I enhancer element of the human interleukin-4 promoter. *Gene* 161: 271-275.

CHROMOSOMAL LOCATION

Genetic locus: CEBPD (human) mapping to 8q11.21.

SOURCE

C/EBP δ (92.69) is a mouse monoclonal antibody raised against recombinant C/EBP δ of human origin.

PRODUCT

Each vial contains 200 μ g IgG₁ kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin. Also available as TransCruz reagent for Gel Supershift and ChIP applications, sc-135733 X, 200 μ g/0.1 ml.

APPLICATIONS

C/EBP δ (92.69) is recommended for detection of C/EBP δ of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)] and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for C/EBP δ siRNA (h): sc-37722, C/EBP δ shRNA Plasmid (h): sc-37722-SH and C/EBP δ shRNA (h) Lentiviral Particles: sc-37722-V.

C/EBP δ (92.69) X TransCruz antibody is recommended for Gel Supershift and ChIP applications.


Molecular Weight of C/EBP δ : 28 kDa.

Positive Controls: A549 cell lysate: sc-2413, HeLa whole cell lysate: sc-2200 or Caco-2 cell lysate: sc-2262.


RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml).

DATA



C/EBP δ (92.69): sc-135733. Western blot analysis of C/EBP δ expression in HeLa (A), Caco-2 (B) and PC-12 (C) whole cell lysates.




C/EBP δ (92.69): sc-135733. Western blot analysis of C/EBP δ expression in A549 whole cell lysate.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures. Not for resale.



See **C/EBP δ (C-6): sc-365546** for C/EBP δ antibody conjugates, including AC, HRP, FITC, PE, Alexa Fluor® 488 and Alexa Fluor® 647.