

Fucokinase (L-12): sc-136692

BACKGROUND

Blood group antigen recognition, metastasis and inflammation utilize fucose, a sugar found in glycoproteins and glycolipids. Fucokinase, also designated L-fucose kinase or FUK, is a 1,084 amino acid enzyme that plays a role in the salvage pathway of fucose reutilization. A member of the GHMP kinase family, Fucokinase exists as two alternative splice variants and catalyzes L-fucose phosphorylation to form β -L-fucose 1-phosphate. The gene encoding Fucokinase maps to human chromosome 16, which encodes over 900 genes and comprises nearly 3% of the human genome. The GAN gene is located on chromosome 16 and, with mutation, may lead to giant axonal neuropathy, a nervous system disorder characterized by increasing malfunction with growth. The rare disorder Rubinstein-Taybi syndrome is also associated with chromosome 16, as is Crohn's disease, which is a gastrointestinal inflammatory condition.

REFERENCES

- Baraitser, M. et al. 1983. The Rubinstein-Taybi syndrome: occurrence in two sets of identical twins. *Clin. Genet.* 23: 318-320.
- Breuning, M.H., et al. 1993. Rubinstein-Taybi syndrome caused by submicroscopic deletions within 16p13.3. *Am. J. Hum. Genet.* 52: 249-254.
- Bomont, P., et al. 2000. The gene encoding gigaxonin, a new member of the cytoskeletal BTB/Kelch repeat family, is mutated in giant axonal neuropathy. *Nat. Genet.* 26: 370-374.
- Hinderlich, S., et al. 2002. Identification of human L-fucose kinase amino acid sequence. *Biochem. Biophys. Res. Commun.* 294: 650-654.
- Kuhlenbäumer, G., et al. 2002. Giant axonal neuropathy (GAN): case report and two novel mutations in the gigaxonin gene. *Neurology* 58: 1273-1276.
- Mathew, C.G. et al. 2004. Genetics of inflammatory bowel disease: progress and prospects. *Hum. Mol. Genet.* 13: R161-R168.
- Online Mendelian Inheritance in Man, OMIM™. 2004. Johns Hopkins University, Baltimore, MD. MIM Number: 608675. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

CHROMOSOMAL LOCATION

Genetic locus: FUK (human) mapping to 16q22.1; Fuk (mouse) mapping to 8 E1.

SOURCE

Fucokinase (L-12) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of Fucokinase of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-136692 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

Fucokinase (L-12) is recommended for detection of Fucokinase isoforms 1 and 2 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).


Fucokinase (L-12) is also recommended for detection of Fucokinase isoforms 1 and 2 in additional species, including bovine.

Suitable for use as control antibody for Fucokinase siRNA (h): sc-93356, Fucokinase siRNA (m): sc-145269, Fucokinase shRNA Plasmid (h): sc-93356-SH, Fucokinase shRNA Plasmid (m): sc-145269-SH, Fucokinase shRNA (h) Lentiviral Particles: sc-93356-V and Fucokinase shRNA (m) Lentiviral Particles: sc-145269-V.

Molecular Weight of Fucokinase: 118 kDa.

Positive Controls: Fucokinase (h): 293T Lysate: sc-115198.

DATA



Fucokinase (L-12): sc-136692. Western blot analysis of Fucokinase expression in non-transfected: sc-117752 (A) and human Fucokinase transfected: sc-115198 (B) 293T whole cell lysates.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.



Try **Fucokinase (F-9): sc-377371**, our highly recommended monoclonal alternative to Fucokinase (L-12).