

# CLLU1OS (E-15): sc-138206

## BACKGROUND

Chronic lymphocytic leukemia (CLL) is an incurable disease characterized by the presence of small mature lymphocytes, intense accumulation of monoclonal B cells and a characteristic CD5 and CD19 co-expression phenotype. While the clinical course of CLL can be highly variable, a CLL specific protein known as CLLU1 (chronic lymphocytic leukemia up-regulated 1) is expressed in CLL patients, with high CLLU1 expression associated with shorter overall survival. Consisting of 121 amino acids, CLLU1 has been predicted to interact with interleukin 4 receptor (IL-4R) and is encoded by a gene located on human chromosome 12. CLLU1OS (chronic lymphocytic leukemia up-regulated 1 opposite strand) is a 101 amino acid protein encoded by a gene that maps to human chromosome 12q22.

## REFERENCES

- Delgado Carrasco, J., Casanova Morcillo, A., Zabalza Alvillos, M. and Ayala Garces, A. 2001. Achondrogenesis type II-hypochondrogenesis: radiological features. *Case report. An. Esp. Pediatr.* 55: 553-557.
- Buhl, A.M., Jurlander, J., Geisler, C.H., Pedersen, L.B., Andersen, M.K., Josefsson, P., Petersen, J.H. and Leffers, H. 2006. CLLU1 expression levels predict time to initiation of therapy and overall survival in chronic lymphocytic leukemia. *Eur. J. Haematol.* 76: 455-464.
- Josefsson, P., Geisler, C.H., Leffers, H., Petersen, J.H., Andersen, M.K., Jurlander, J. and Buhl, A.M. 2007. CLLU1 expression analysis adds prognostic information to risk prediction in chronic lymphocytic leukemia. *Blood* 109: 4973-4979.
- Chen, L., Li, J., Zheng, W., Zhang, Y., Wu, Y., Li, L., Qian, S. and Xu, W. 2007. The prognostic evaluation of CLLU1 expression levels in 50 Chinese patients with chronic lymphocytic leukemia. *Leuk. Lymphoma* 48: 1785-1792.
- Chen, L.J., Zheng, W.J., Wu, Y.J., Li, L., Qian, S.X., Xu, W. and Li, J.Y. 2007. Expression of CLLU1 in patients with chronic lymphocytic leukemia and its prognostic significance. *Zhonghua Xue Ye Xue Za Zhi* 28: 737-740.
- Lo, F.S., Luo, J.D., Lee, Y.J., Shu, S.G., Kuo, M.T. and Chiou, C.C. 2009. High resolution melting analysis for mutation detection for PTPN11 gene: applications of this method for diagnosis of Noonan syndrome. *Clin. Chim. Acta* 409: 75-77.
- Benussi, D.G., Costa, P., Zollino, M., Murdolo, M., Petix, V., Carrozzi, M. and Pecile, V. 2009. Trisomy 12p and monosomy 4p: phenotype-genotype correlation. *Genet. Test. Mol. Biomarkers* 13: 199-204.
- Buhl, A.M., Novotny, G.W., Josefsson, P., Nielsen, J.E., Pedersen, L.B., Geisler, C., Rassenti, L.Z., Kipps, T.J., Jurlander, J. and Leffers, H. 2009. The CLLU1 expression level is a stable and inherent feature of the chronic lymphocytic leukemia clone. *Leukemia* 23: 1182-1186.

## CHROMOSOMAL LOCATION

Genetic locus: CLLU1OS (human) mapping to 12q22.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## SOURCE

CLLU1OS (E-15) is an affinity purified rabbit polyclonal antibody raised against a peptide mapping within an internal region of CLLU1OS of human origin.

## PRODUCT

Each vial contains 100 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-138206 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

CLLU1OS (E-15) is recommended for detection of CLLU1OS of human origin by Western Blotting (starting dilution 1:100, dilution range 1:50-1:500), immunofluorescence (starting dilution 1:25, dilution range 1:25-1:250) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with CLLU1.

Suitable for use as control antibody for CLLU1OS siRNA (h): sc-96183, CLLU1OS shRNA Plasmid (h): sc-96183-SH and CLLU1OS shRNA (h) Lentiviral Particles: sc-96183-V.

Molecular Weight of CLLU1OS: 11 kDa.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.