

SPATA5 (C-12): sc-138638

BACKGROUND

SPATA5 (spermatogenesis-associated protein 5), also known as SPAF (spermatogenesis-associated factor protein) or AFG2 (ATPase family protein 2 homolog), is an 893 amino acid protein that localizes to cytoplasm and mitochondrion, and may be involved in morphological and functional mitochondrial transformations during spermatogenesis. Existing as three alternatively spliced isoforms, SPATA5 belongs to the AAA ATPase family and the AFG2 subfamily. The gene that encodes SPATA5 consists of more than 396,000 bases and maps to human chromosome 4q28.1. Housing nearly 900 genes, chromosome 4 represents approximately 6% of the human genome and is associated with Huntington's disease, Ellis-van Creveld syndrome, methylmalonic acidemia and polycystic kidney disease.

REFERENCES

1. Bonaventure, J., et al. 1996. Common mutations in the fibroblast growth factor receptor 3 (FGFR 3) gene account for achondroplasia, hypochondroplasia, and thanatophoric dwarfism. *Am. J. Med. Genet.* 63: 148-154.
2. Kalchman, M.A., et al. 1996. Huntingtin is ubiquitinated and interacts with a specific ubiquitin-conjugating enzyme. *J. Biol. Chem.* 271: 19385-19394.
3. Howard, T.D., et al. 1997. Autosomal dominant postaxial polydactyly, nail dystrophy, and dental abnormalities map to chromosome 4p16, in the region containing the Ellis-van Creveld syndrome locus. *Am. J. Hum. Genet.* 61: 1405-1412.
4. Singhrao, S.K., et al. 1998. Huntingtin protein colocalizes with lesions of neurodegenerative diseases: An investigation in Huntington's, Alzheimer's, and Pick's diseases. *Exp. Neurol.* 150: 213-222.
5. Krakow, D., et al. 2000. Exclusion of the Ellis-van Creveld region on chromosome 4p16 in some families with asphyxiating thoracic dystrophy and short-rib polydactyly syndromes. *Eur. J. Hum. Genet.* 8: 645-648.
6. Liu, Y., et al. 2000. SPAF, a new AAA-protein specific to early spermatogenesis and malignant conversion. *Oncogene* 19: 1579-1588.
7. Sommardahl, C., et al. 2001. Phenotypic variations of orpk mutation and chromosomal localization of modifiers influencing kidney phenotype. *Physiol. Genomics* 7: 127-134.
8. Dobson, C.M., et al. 2002. Identification of the gene responsible for the cblA complementation group of vitamin B12-responsive methylmalonic acidemia based on analysis of prokaryotic gene arrangements. *Proc. Natl. Acad. Sci. USA* 99: 15554-15559.
9. Online Mendelian Inheritance in Man, OMIM[™]. 2011. Johns Hopkins University, Baltimore, MD. MIM Number: 613940. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

CHROMOSOMAL LOCATION

Genetic locus: SPATA5 (human) mapping to 4q28.1; Spata5 (mouse) mapping to 3 B.

SOURCE

SPATA5 (C-12) is an affinity purified rabbit polyclonal antibody raised against a peptide mapping at the C-terminus of SPATA5 of human origin.

PRODUCT

Each vial contains 100 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-138638 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

SPATA5 (C-12) is recommended for detection of SPATA5 isoform 1 of mouse, rat and human origin by Western Blotting (starting dilution 1:100, dilution range 1:50-1:500), immunofluorescence (starting dilution 1:25, dilution range 1:25-1:250) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other SPATA family members or isoforms SPATA5-2 or SPATA5-3.

SPATA5 (C-12) is also recommended for detection of SPATA5 isoform 1 in additional species, including equine, bovine and porcine.

Suitable for use as control antibody for SPATA5 siRNA (h): sc-89261, SPATA5 siRNA (m): sc-153722, SPATA5 shRNA Plasmid (h): sc-89261-SH, SPATA5 shRNA Plasmid (m): sc-153722-SH, SPATA5 shRNA (h) Lentiviral Particles: sc-89261-V and SPATA5 shRNA (m) Lentiviral Particles: sc-153722-V.

Molecular Weight of SPATA5 isoforms: 98/86/76 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker[™] compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker[™] Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz[™] Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.