

WDR78 (V-13): sc-138732

BACKGROUND

WD-repeats are motifs that are found in a variety of proteins and are characterized by a conserved core of 40-60 amino acids that commonly form a tertiary propeller structure. While proteins that contain WD-repeats participate in a wide range of cellular functions, they are generally involved in regulatory mechanisms concerning chromatin assembly, cell cycle control, signal transduction, RNA processing, apoptosis and vesicular trafficking. WDR78 (WD repeat domain 78) is an 848 amino acid protein that contains 6 WD repeats and exists as three alternatively spliced isoforms. The gene encoding WDR78 maps to human chromosome 1, which spans 260 million base pairs, contains over 3,000 genes and comprises nearly 8% of the human genome.

REFERENCES

1. van der Voorn, L. and Ploegh, H.L. 1992. The WD-40 repeat. *FEBS Lett.* 307: 131-134.
2. Neer, E.J., Schmidt, C.J., Nambudripad, R. and Smith, T.F. 1994. The ancient regulatory-protein family of WD-repeat proteins. *Nature* 371: 297-300.

CHROMOSOMAL LOCATION

Genetic locus: WDR78 (human) mapping to 1p31.3; Wdr78 (mouse) mapping to 4 C6.

SOURCE

WDR78 (V-13) is an affinity purified rabbit polyclonal antibody raised against a peptide mapping within an internal region of WDR78 of human origin.

PRODUCT

Each vial contains 100 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-138732 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

WDR78 (V-13) is recommended for detection of WDR78 isoforms 1 and 3 of mouse, rat and human origin by Western Blotting (starting dilution 1:100, dilution range 1:50-1:500), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:25, dilution range 1:25-1:250) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other WDR family members or WDR78 isoform 2.

WDR78 (V-13) is also recommended for detection of WDR78 isoforms 1 and 3 in additional species, including bovine.

Suitable for use as control antibody for WDR78 siRNA (h): sc-78584, WDR78 siRNA (m): sc-155318, WDR78 shRNA Plasmid (h): sc-78584-SH, WDR78 shRNA Plasmid (m): sc-155318-SH, WDR78 shRNA (h) Lentiviral Particles: sc-78584-V and WDR78 shRNA (m) Lentiviral Particles: sc-155318-V.


Molecular Weight of WDR78 isoforms: 95/57/63 kDa.

Positive Controls: WDR78 (m): 293T Lysate: sc-126238.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



WDR78 (V-13): sc-138732. Western blot analysis of WDR78 expression in non-transfected: sc-117752 (A) and mouse WDR78 transfected: sc-126238 (B) 293T whole cell lysates.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

MONOS
Satisfaction
Guaranteed

Try **WDR78 (E-2): sc-390401** or **WDR78 (E-7): sc-390633**, our highly recommended monoclonal alternatives to WDR78 (V-13).