

HB-EGF (M-18): sc-1414

BACKGROUND

Heparin binding epidermal-like growth factor (HB-EGF), a member of the EGF family of mitogens, binds to the EGF receptor (EGFR) and to heparin sulfate proteoglycans on the cell surface. HB-EGF was originally isolated from medium conditioned by the growth of the human histocytic lymphoma cell U-937 on the basis of its heparin-binding ability and its mitogenic activity for Balb-3T3 fibroblasts. The HB-EGF gene encodes a 208 amino acid precursor containing a signal peptide and transmembrane domain. Mature HB-EGF is a soluble protein 86 amino acids in length and results from the enzymatic cleavage of the membrane bound precursor. The membrane-bound form of HB-EGF has been identified as the diphtheria toxin receptor. Preincubation of Vero cells with phorbol 12-myristate 13-acetate (PMA) induces the proteolytic cleavage of HB-EGF outside the membrane anchor.

CHROMOSOMAL LOCATION

Genetic locus: HBEGF (human) mapping to 5q31.3; Hbegf (mouse) mapping to 18 B2.

SOURCE

HB-EGF (M-18) is an affinity purified goat polyclonal antibody raised against a peptide mapping at a C-terminal cytoplasmic domain of HB-EGF of mouse origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-1414 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

HB-EGF (M-18) is recommended for detection of precursor HB-EGF of mouse, rat and, to a lesser extent, human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for HB-EGF siRNA (h): sc-39420, HB-EGF siRNA (m): sc-39421, HB-EGF shRNA Plasmid (h): sc-39420-SH, HB-EGF shRNA Plasmid (m): sc-39421-SH, HB-EGF shRNA (h) Lentiviral Particles: sc-39420-V and HB-EGF shRNA (m) Lentiviral Particles: sc-39421-V.

Molecular Weight of HB-EGF: 22 kDa.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.


PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

DATA



HB-EGF (M-18): sc-1414. Western blot analysis of human recombinant HB-EGF.

SELECT PRODUCT CITATIONS

- Hamilton, M. and Wolfman, A. 1998. Oncogenic Ha-Ras-dependent mitogen-activated protein kinase activity requires signaling through the epidermal growth factor receptor. *J. Biol. Chem.* 273: 28155-28162.
- Ku, S.Y., et al. 2005. Gonadotropin stimulation may not affect the expression of heparin-binding epidermal growth factor-like growth factor in the murine endometrium. *Gynecol. Obstet. Invest.* 59: 179-183.
- Kaya, G., et al. 2006. Hyaluronate fragments reverse skin atrophy by a CD44-dependent mechanism. *PLoS Med.* 3: e493.
- Chen, H., et al. 2007. Epidermal growth factor receptor in adult retinal neurons of rat, mouse, and human. *J. Comp. Neurol.* 500: 299-310.
- Charalambous, C.T., et al. 2007. Latent membrane protein 1-induced EGFR signalling is negatively regulated by TGF α prior to neoplasia. *Carcinogenesis* 28: 1839-1848.
- Ortiz, C., et al. 2008. Inhibition of the EGF receptor blocks autocrine growth and increases the cytotoxic effects of doxorubicin in rat hepatoma cells: role of reactive oxygen species production and glutathione depletion. *Biochem. Pharmacol.* 75: 1935-1945.
- Torres, C., et al. 2008. A cyclic peptide derived from α -fetoprotein inhibits the proliferative effects of the epidermal growth factor and estradiol in MCF7 cells. *Oncol. Rep.* 19: 1597-1603.
- Edwards, J.P., et al. 2009. The expression of heparin-binding epidermal growth factor-like growth factor by regulatory macrophages. *J. Immunol.* 182: 1929-1939.
- Barnes, L., et al. 2010. Synergistic effect of hyaluronate fragments in retinaldehyde-induced skin hyperplasia which is a Cd44-dependent phenomenon. *PLoS ONE* 5: e14372.
- Hamaoka, M., et al. 2010. Anti-human HB-EGF monoclonal antibodies inhibiting ectodomain shedding of HB-EGF and diphtheria toxin binding. *J. Biochem.* 148: 55-69.
- Ozcakir, H.T., et al. 2011. Effect of GnRH antagonist therapy on the expression of MUC-1 and heparin binding growth factor expression in the endometrium of hyperstimulated rats. *Clin. Exp. Obstet. Gynecol.* 38: 76-80.