

FANCE siRNA (m): sc-145063

BACKGROUND

Fanconi anemia (FA) is an autosomal recessive disorder characterized by bone marrow failure, birth defects and chromosomal instability. At the cellular level, FA is characterized by spontaneous chromosomal breakage and a unique hypersensitivity to DNA cross-linking agents. At least eight complementation groups (A-G) have been identified and six FA genes (for subtypes A, C, D2, E, F, and G) have been cloned. The FA proteins lack sequence homologies or motifs that could point to a molecular function. Phosphorylation of FANCF (Fanconi anemia complementation group) proteins is thought to be important for the function of the FA pathway. FA proteins encoded by six cloned FA genes (FANCA, FANCC, FANCD2, FANCE, FANCF and FANCG) cooperate in a common pathway, culminating in the monoubiquitination of FANCD2 protein and the co-localization of FANCD2 and BRCA1 proteins in nuclear foci. The human FANCE gene maps to chromosome 6p21.31, contains 10 exons and encodes a novel 536 amino acid protein with two potential nuclear localization signals.

REFERENCES

1. de Winter, J.P., et al. 2000. The Fanconi anemia protein FANCF forms a nuclear complex with FANCA, FANCC and FANCG. *Hum. Mol. Genet.* 9: 2665-2674.
2. de Winter, J.P., et al. 2000. Isolation of a cDNA representing the Fanconi anemia complementation group E gene. *Am. J. Hum. Genet.* 67: 1306-1308.
3. Yagasaki, H., et al. 2001. A cytoplasmic serine protein kinase binds and may regulate the Fanconi anemia protein FANCA. *Blood* 98: 3650-3657.
4. Wilson, J.B., et al. 2001. The Chinese hamster FANCG/XRCC9 mutant NM3 fails to express the monoubiquitinated form of the FANCD2 protein, is hypersensitive to a range of DNA damaging agents and exhibits a normal level of spontaneous sister chromatid exchange. *Carcinogenesis* 22: 1939-1946.
5. Siddique, M.A., et al. 2001. Function of the Fanconi anemia pathway in Fanconi anemia complementation group F and D1 cells. *Exp. Hematol.* 29: 1448-1455.
6. Otsuki, T., et al. 2002. Phosphorylation of Fanconi anemia protein, FANCA, is regulated by Akt kinase. *Biochem. Biophys. Res. Commun.* 291: 628-634.
7. Taniguchi, T. and Dandrea, A.D. 2002. Molecular pathogenesis of fanconi anemia. *Int. J. Hematol.* 75: 123-128.

CHROMOSOMAL LOCATION

Genetic locus: Fance (mouse) mapping to 17 A3.3.

PRODUCT

FANCE siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see FANCE shRNA Plasmid (m): sc-145063-SH and FANCE shRNA (m) Lentiviral Particles: sc-145063-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

FANCE siRNA (m) is recommended for the inhibition of FANCE expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

FANCE (B-1): sc-398558 is recommended as a control antibody for monitoring of FANCE gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker[™] Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz[®] Mounting Medium: sc-24941 or UltraCruz[®] Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor FANCE gene expression knockdown using RT-PCR Primer: FANCE (m)-PR: sc-145063-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.