

glypican-1 (C-18): sc-14648

BACKGROUND

Glypican-1 (GPC1), a member of the glycosylphosphatidylinositol anchored cell surface heparan sulfate proteoglycans, is involved with cell adhesion and migration, lipoprotein metabolism, modulation of growth factor activities and anticoagulation. Glypican-1 binds to and modulates the activity of several fibroblast growth factors (FGFs) including FGF-1, FGF-2 and FGF-7. Glypican-1 acts as an extracellular chaperone for VEGF165 to help restore receptor binding ability after oxidation. The heparan sulfate chains of glypican-1 mediate specific binding of glypican-1 to VEGF165. When present on the surface of marrow stromal cells, glypican-1 may aid in the maintenance and development of hematopoietic stem and progenitor cells. Human pancreatic cancer cells express a large amount of glypican-1 when compared to glypican-1 levels in normal pancreatic cells. Glypican-1 may play an important role in the response of pancreatic cancer cells to mitogenic stimuli, such as FGF-2. The gene encoding human glypican-1 maps to chromosome 2q37.3.

REFERENCES

- David, G. 1993. Integral membrane heparan sulfate proteoglycans. *FASEB J.* 7: 1023-1030.
- Vermeesch, J.R., et al. 1995. Assignment of the human glypican gene (GPC1) to 2q35-q37 by fluorescence *in situ* hybridization. *Genomics* 25: 327-339.
- Steinfeld, R., et al. 1996. Stimulation of fibroblast growth factor receptor-1 occupancy and signaling by cell surface-associated syndecans and glypican. *J. Cell Biol.* 133: 405-416.
- Weksberg, R., et al. 1996. Glypicans: a growing trend. *Nat. Genet.* 12: 225-227.
- Bonneh-Barkay, D., et al. 1997. Identification of glypican as a dual modulator of the biological activity of fibroblast growth factors. *J. Biol. Chem.* 272: 12415-12421.
- Kleeff, J., et al. 1998. The cell-surface heparan sulfate proteoglycan glypican-1 regulates growth factor action in pancreatic carcinoma cells and is overexpressed in human pancreatic cancer. *J. Clin. Invest.* 102: 1662-1673.
- Gengrinovitch, S., et al. 1999. Glypican-1 is a VEGF165 binding proteoglycan that acts as an extracellular chaperone for VEGF165. *J. Biol. Chem.* 274: 10816-10822.
- Schofield, K.P., et al. 1999. Expression of proteoglycan core proteins in human bone marrow stroma. *Biochem. J.* 3: 663-668.

CHROMOSOMAL LOCATION

Genetic locus: GPC1 (human) mapping to 2q37.3; Gpc1 (mouse) mapping to 1 D.

SOURCE

glypican-1 (C-18) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of glypican-1 of human origin.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-14648 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

glypican-1 (C-18) is recommended for detection of glypican-1 of mouse, rat and human Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

glypican-1 (C-18) is also recommended for detection of glypican-1 in additional species, including bovine and porcine.

Suitable for use as control antibody for glypican-1 siRNA (h): sc-40638, glypican-1 siRNA (m): sc-40639, glypican-1 shRNA Plasmid (h): sc-40638-SH, glypican-1 shRNA Plasmid (m): sc-40639-SH, glypican-1 shRNA (h) Lentiviral Particles: sc-40638-V and glypican-1 shRNA (m) Lentiviral Particles: sc-40639-V.

Molecular Weight of glypican-1: 64 kDa.

Positive Controls: K-562 whole cell lysate: sc-2203, mouse brain extract: sc-2253 or LADMAC whole cell lysate: sc-364189.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.



Try **glypican-1 (A-10): sc-365000** or **glypican-1 (4D1): sc-101827**, our highly recommended monoclonal alternatives to glypican-1 (C-18).