

EYA4 (E-20): sc-15107

BACKGROUND

A gene of chromosome 6q23.2 encodes the 640 amino acid protein, EYA4 (eyes absent). EYA is one of four members of the eyes absent family. A 271 amino acid domain at the carboxyl terminal is highly conserved amongst the members of the eyes absent family. EYA4 is expressed in the craniofacial mesenchyme, the dermamyotome, and the limb. The conserved region in other EYA proteins interacts with SIX, DACH, and G-proteins, which regulate transcription in early embryonic development. SIX translocates EYA1-3 to the nucleus, and G proteins can stop this interaction. Premature stop codon mutations in EYA4 cause postlingual, progressive autosomal dominant hearing loss in humans. This shows that EYA4 is also vital to the mature organ of Corti. EYA4 may cause oculo-dento-digital syndrome, based on its expression pattern and map position.

REFERENCES

1. Borsani, G., DeGrandi, A., Ballabio, A., Bulfone, A., Bernard, L., Banfi, S., Gattuso, C., Mariani, M., Dixon, M., Donnai, D., Metcalfe, K., Winter, R., Robertson, M., Axton, R., Brown, A., van Heyningen, V. and Hanson, I. 1999. EYA4, a novel vertebrate gene related to *Drosophila* eyes absent. *Hum. Molec. Genet.* 8: 11-23.
2. Ohto, H., Kamada, S., Tago, K., Tominaga, S.I., Ozaki, H., Sato, S. and Kawakami, K. 1999. Cooperation of six and eya in activation of their target genes through nuclear translocation of Eya. *Mol. Cell. Biol.* 19: 6815-6824.
3. Heanue, T.A., Reshef, R., Davis, R.J., Mardon, G., Oliver, G., Tomarev, S., Lassar, A.B. and Tabin, C.J. 1999. Synergistic regulation of vertebrate muscle development by Dach2, Eya2, and Six1, homologs of genes required for *Drosophila* eye formation. *Genes Dev.* 13: 3231-3243.
4. Fan, X., Brass, L.F., Poncz, M., Spitz, F., Maire, P. and Manning, D.R. 2000. The α subunits of G₂ and G_i interact with the eyes absent transcription cofactor Eya2, preventing its interaction with the six class of homeodomain-containing proteins. *J. Biol. Chem.* 275: 32129-32134.
5. Wayne, S., Robertson, N.G., DeClau, F., Chen, N., Verhoeven, K., Prasad, S., Tranebjarg, L., Morton, C.C., Ryan, A.F., Van Camp, G. and Smith, R.J. 2001. Mutations in the transcriptional activator EYA4 cause late-onset deafness at the DFNA10 locus. *Hum. Mol. Genet.* 10: 195-200.

CHROMOSOMAL LOCATION

Genetic locus: EYA4 (human) mapping to 6q23.2; Eya4 (mouse) mapping to 10 A3.

SOURCE

EYA4 (E-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of EYA4 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-15106 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

EYA4 (E-20) is recommended for detection of EYA4 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

EYA4 (E-20) is also recommended for detection of EYA4 in additional species, including canine, bovine, porcine and avian.

Suitable for use as control antibody for EYA4 siRNA (h): sc-41952, EYA4 siRNA (m): sc-41952, EYA4 shRNA Plasmid (h): sc-41952-SH, EYA4 shRNA Plasmid (m): sc-41952-SH, EYA4 shRNA (h) Lentiviral Particles: sc-41952-V and EYA4 shRNA (m) Lentiviral Particles: sc-41952-V.

Molecular Weight (predicted) of EYA4: 70 kDa.

Molecular Weight (observed) of EYA4: 86-90 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

MONOS
Satisfaction
Guaranteed

Try **EYA4 (E-11): sc-393111**, our highly recommended monoclonal alternative to EYA4 (E-20).