

TPH (D-15): sc-15114

BACKGROUND

Phenylalanine hydroxylase (PAH), tyrosine hydroxylase (TH) and tryptophan hydroxylase (TPH) comprise a small family of monooxygenases that use tetrahydropterine as a cofactor during the catabolism of aromatic L-amino acids. PAH, TH and TPH all contain catalytic domains with an amino-terminal regulatory domain and a short carboxy-terminal tetramerization domain. Each of these enzymes also contains a single ferrous iron atom, which is bound to two histidines and a glutamate, and is likely to be involved in the formation of the hydroxylating intermediate. TPH is both the first and rate-limiting-step in the biosynthesis of serotonin in the central nervous system and melatonin in the pineal gland. Alteration of TPH function may be a key factor in the pathology of several neuropsychiatric disorders associated with serotonin, including depression, aggression, alcoholism and schizophrenia. For instance, L-DOPA, which is used as a common therapy for Parkinson's disease (PD) patients, inhibits TPH function which, subsequently, is thought to contribute to the onset of depression in PD patients.

REFERENCES

1. Abbar, M., Courtet, P., Malafosse, A. and Castelnau, D. 1996. Epidemiologic and molecular genetic of suicidal behavior. *Encephale* 22: 19-24.
2. Mockus, S.M. and Vrana, K.E. 1998. Advances in the molecular characterization of tryptophan hydroxylase. *J. Mol. Neurosci.* 10: 163-179.
3. Fitzpatrick, P.F. 1999. Tetrahydropterin-dependent amino acid hydroxylases. *Annu. Rev. Biochem.* 68: 355-381.
4. Kuhn, D.M. 1999. Tryptophan hydroxylase regulation. Drug-induced modifications that alter serotonin neuronal function. *Adv. Exp. Med. Biol.* 467: 19-27.
5. Kowlessur, D. and Kaufman, S. 1999. Cloning and expression of recombinant human pineal tryptophan hydroxylase in *Escherichia coli*: purification and characterization of the cloned enzyme. *Biochim. Biophys. Acta* 1434: 317-330.
6. Nagatsu, T. and Ichinose, H. 1999. Regulation of pteridine-requiring enzymes by the cofactor tetrahydrobiopterin. *Mol. Neurobiol.* 19: 79-96.
7. Veenstra-VanderWeele, J., Anderson, G.M. and Cook, E.H.J. 2000. Pharmacogenetics and the serotonin system: initial studies and future directions. *Eur. J. Pharmacol.* 410: 165-181.

CHROMOSOMAL LOCATION

Genetic locus: TPH1 (human) mapping to 11p15.1; Tph1 (mouse) mapping to 7 B4.

SOURCE

TPH (D-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of TPH of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-15114 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

TPH (D-15) is recommended for detection of TPH of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

PH (D-15) is also recommended for detection of TPH in additional species, including equine and canine.

Suitable for use as control antibody for TPH siRNA (h): sc-41526, TPH siRNA (m): sc-41527, TPH shRNA Plasmid (h): sc-41526-SH, TPH shRNA Plasmid (m): sc-41527-SH, TPH shRNA (h) Lentiviral Particles: sc-41526-V and TPH shRNA (m) Lentiviral Particles: sc-41527-V.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.