

# GDNF siRNA (r): sc-156116

## BACKGROUND

Glial cell line-derived neurotrophic factor (GDNF) has been identified as a potent neurotrophic factor that enhances survival of midbrain dopaminergic neurons. GDNF is a glycosylated, disulfide-bonded homodimer and is a distantly related member of the TGF $\beta$  superfamily of growth regulatory ligands. GDNF contains the seven conserved cysteine residues in the same relative spacing characteristic of all members of the TGF $\beta$  superfamily. In embryonic midbrain cultures, GDNF promotes the survival and morphological differentiation of dopaminergic neurons and increases their high-affinity dopamine uptake. On the basis of these findings, it has been suggested that GDNF may have utility in the treatment of Parkinson's disease, which is marked by progressive degeneration of midbrain dopaminergic neurons.

## REFERENCES

- Schubert, D., et al. 1974. Clonal cell lines from the rat central nervous system. *Nature* 249: 224-227.
- Derynck, R., et al. 1985. Human transforming growth factor  $\beta$  cDNA sequence and expression in tumor cell lines. *Nature* 316: 701-705.
- ten Dijke, P., et al. 1988. Identification of a new member of the transforming growth factor type  $\beta$  gene family. *Proc. Natl. Acad. Sci. USA* 85: 4715-4719.
- Miller, D.A., et al. 1990. Transforming growth factor  $\beta$ : a family of growth regulatory peptides. *Ann. N.Y. Acad. Sci.* 593: 208-217.
- Lin, L.H., et al. 1993. GDNF: a glial cell line-derived neurotrophic factor for midbrain dopaminergic neurons. *Science* 260: 1130-1132.
- Pozas, E., et al. 2005. GDNF and GFR $\alpha$ -1 promote differentiation and tangential migration of cortical GABAergic neurons. *Neuron* 45: 701-713.

## CHROMOSOMAL LOCATION

Genetic locus: *Gdnf* (rat) mapping to 2q16.

## PRODUCT

GDNF siRNA (r) is a pool of 2 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see GDNF shRNA Plasmid (r): sc-156116-SH and GDNF shRNA (r) Lentiviral Particles: sc-156116-V as alternate gene silencing products.

For independent verification of GDNF (r) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-156116A and sc-156116B.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20 $^{\circ}$  C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20 $^{\circ}$  C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

GDNF siRNA (r) is recommended for the inhibition of GDNF expression in rat cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## GENE EXPRESSION MONITORING

GDNF (B-8): sc-13147 is recommended as a control antibody for monitoring of GDNF gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG $\kappa$  BP-HRP: sc-516102 or m-IgG $\kappa$  BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker<sup>TM</sup> Molecular Weight Standards: sc-2035, UltraCruz<sup>®</sup> Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG $\kappa$  BP-FITC: sc-516140 or m-IgG $\kappa$  BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz<sup>®</sup> Mounting Medium: sc-24941 or UltraCruz<sup>®</sup> Hard-set Mounting Medium: sc-359850.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor GDNF gene expression knockdown using RT-PCR Primer: GDNF (r)-PR: sc-156116-PR (20  $\mu$ l, 460 bp). Annealing temperature for the primers should be 55-60 $^{\circ}$  C and the extension temperature should be 68-72 $^{\circ}$  C.

## SELECT PRODUCT CITATIONS

- Han, T.Y., et al. 2015. Intestinal smooth muscle phenotype determines enteric neuronal survival via GDNF expression. *Neuroscience* 290: 357-368.
- Lv, B., et al. 2016. Activated microglia induce bone marrow mesenchymal stem cells to produce glial cell-derived neurotrophic factor and protect neurons against oxygen-glucose deprivation injury. *Front. Cell. Neurosci.* 10: 283.
- Zhong, Z., et al. 2020. Bone marrow mesenchymal stem cells upregulate PI3K/Akt pathway and down-regulate NF $\kappa$ B pathway by secreting glial cell-derived neurotrophic factors to regulate microglial polarization and alleviate deafferentation pain in rats. *Neurobiol. Dis.* E-published.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.