

T3JAM (S-12): sc-160101

BACKGROUND

T3JAM (TRAF3 interacting protein 3), also known as TRAF3-interacting JNK-activating modulator or TRAF3IP3, is a 551 amino acid single-pass type IV membrane protein that acts as an adaptor molecule and stimulates cell growth by modulating the c-Jun N-terminal kinase (JNK) pathway in conjunction with TRAF3. Expressed in thymus, bone marrow and spleen, T3JAM binds the isoleucine zipper of TRAF3 with its coiled-coil domain. T3JAM has been found to associate specifically with TRAF3 but not with other TRAF family members. As a result of alternative splicing events, three T3JAM isoforms exist. The gene encoding T3JAM maps to human chromosome 1, which spans about 260 million base pairs and comprises nearly 8% of the human genome. Aberrations in chromosome 1 are found in a variety of cancers including head and neck cancer, malignant melanoma and multiple myeloma.

REFERENCES

1. Dadgostar, H., Doyle, S.E., Shahangian, A., Garcia, D.E. and Cheng, G. 2003. T3JAM, a novel protein that specifically interacts with TRAF3 and promotes the activation of JNK1. *FEBS Lett.* 553: 403-407.
2. Zhang, C., Xie, L., Cheng, H. and Wang, Y. 2007. TRAF3 interacts with Smac/DIABLO and enhances the proapoptotic effect of Smac/DIABLO in cytoplasm. *Acta Biochim. Biophys. Sin.* 39: 108-116.
3. Ma, X., Wang, X., Gao, X., Wang, L., Lu, Y., Gao, P., Deng, W., Yu, P., Ma, J., Guo, J., Cheng, H., Zhang, C., Shi, T. and Ma, D. 2007. Identification of five human novel genes associated with cell proliferation by cell-based screening from an expressed cDNA ORF library. *Life Sci.* 81: 1141-1151.
4. Balcárková, J., Urbánková, H., Scudla, V., Holzerová, M., Bacovský, J., Indrák, K. and Jarosová, M. 2009. Gain of chromosome arm 1q in patients in relapse and progression of multiple myeloma. *Cancer Genet. Cytogenet.* 192: 68-72.
5. Online Mendelian Inheritance in Man, OMIM™. 2007. Johns Hopkins University, Baltimore, MD. MIM Number: 608255. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

CHROMOSOMAL LOCATION

Genetic locus: TRAF3IP3 (human) mapping to 1q32.2; Traf3ip3 (mouse) mapping to 1 H6.

SOURCE

T3JAM (S-12) is an affinity purified goat polyclonal antibody raised against a peptide mapping within a cytoplasmic domain of T3JAM of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-160101 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

RESEARCH USE

For research use only, not for use in diagnostic procedures.

APPLICATIONS

T3JAM (S-12) is recommended for detection of T3JAM isoforms 1-3 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for T3JAM siRNA (h): sc-78894, T3JAM siRNA (m): sc-154028, T3JAM shRNA Plasmid (h): sc-78894-SH, T3JAM shRNA Plasmid (m): sc-154028-SH, T3JAM shRNA (h) Lentiviral Particles: sc-78894-V and T3JAM shRNA (m) Lentiviral Particles: sc-154028-V.

Molecular Weight of T3JAM: 64 kDa.

Positive Controls: HEL 92.1.7 cell lysate: sc-2270.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.


 MONOS
Satisfation
Guaranteed

Try **T3JAM (B-7): sc-398781** or **T3JAM (H-2): sc-398895**, our highly recommended monoclonal alternatives to T3JAM (S-12).