

MVD (V-19): sc-160552

BACKGROUND

MVD (mevalonate [diphospho] decarboxylase), also known as MPD (mevalonate pyrophosphate decarboxylase), is a 400 amino acid protein that belongs to the diphosphomevalonate decarboxylase family. Expressed in lung, liver, heart, skeletal muscle, brain, pancreas, placenta and kidney, MVD enzymatically catalyzes the first step in isoprene biosynthesis, namely the ATP-dependent conversion of mevalonate pyrophosphate into isopentenyl pyrophosphate, a cholesterol precursor. MVD exists as a homodimer that simultaneously dehydrates and decarboxylates its substrate while hydrolyzing ATP. As MVD is a crucial enzyme in early cholesterol synthesis, it may be a useful target for drugs that aim to lower serum cholesterol levels.

REFERENCES

1. Ku, E.C. 1996. Regulation of fatty acid biosynthesis by intermediates of the cholesterol biosynthetic pathway. *Biochem. Biophys. Res. Commun.* 225: 173-179.
2. Toth, M.J. and Huwyler, L. 1996. Molecular cloning and expression of the cDNAs encoding human and yeast mevalonate pyrophosphate decarboxylase. *J. Biol. Chem.* 271: 7895-7898.
3. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 603236. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
4. Michihara, A., Akasaki, K., Yamori, Y. and Tsuji, H. 2003. Subcellular distribution of mouse mevalonate pyrophosphate decarboxylase. *Biol. Pharm. Bull.* 26: 579-584.
5. Wadhwa, R., Yaguchi, T., Hasan, M.K., Taira, K. and Kaul, S.C. 2003. Mortalin-MPD (mevalonate pyrophosphate decarboxylase) interactions and their role in control of cellular proliferation. *Biochem. Biophys. Res. Commun.* 302: 735-742.
6. Flock, G., Baggio, L.L., Longuet, C. and Drucker, D.J. 2007. Incretin receptors for glucagon-like peptide 1 and glucose-dependent Insulinotropic polypeptide are essential for the sustained metabolic actions of vildagliptin in mice. *Diabetes* 56: 3006-3013.

CHROMOSOMAL LOCATION

Genetic locus: MVD (human) mapping to 16q24.3; Mvd (mouse) mapping to 8 E1.

SOURCE

MVD (V-19) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of MVD of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-160552 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

MVD (V-19) is recommended for detection of MVD of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

MVD (V-19) is also recommended for detection of MVD in additional species, including equine, canine and bovine.

Suitable for use as control antibody for MVD siRNA (h): sc-93276, MVD siRNA (m): sc-149724, MVD shRNA Plasmid (h): sc-93276-SH, MVD shRNA Plasmid (m): sc-149724-SH, MVD shRNA (h) Lentiviral Particles: sc-93276-V and MVD shRNA (m) Lentiviral Particles: sc-149724-V.

Molecular Weight of MVD: 43 kDa.

Positive Controls: A-431 whole cell lysate: sc-2201 or K-562 whole cell lysate: sc-2203.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.



Try **MVD (H-11): sc-376975** or **MVD (2B5): sc-100559**, our highly recommended monoclonal alternatives to MVD (V-19).