

# NUBP1 (M-13): sc-160614

## BACKGROUND

NUBP1 (nucleotide binding protein 1), also known as cytosolic Fe-S cluster assembly factor NUBP1 or Nbp35, is a 320 amino acid protein involved in the regulation of centrosome duplication. A member of the Mrp/NBP35 ATP-binding protein family and the NUBP1/NBP35 subfamily, NUBP1 exists as two alternatively spliced isoforms and is known to interact with KIFC1 and NUBP2. NUBP1 is a component of the cytosolic iron-sulfur (Fe/S) protein assembly machinery and can transfer iron-sulfur clusters to certain apoproteins. The gene encoding NUBP1 maps to human chromosome 16, which encodes over 900 genes and comprises nearly 3% of the human genome. The GAN gene is located on chromosome 16 and, with mutation, may lead to giant axonal neuropathy, a nervous system disorder characterized by increasing malfunction with growth. The rare disorder Rubinstein-Taybi syndrome is also associated with chromosome 16, as is Crohn's disease, which is a gastrointestinal inflammatory condition.

## REFERENCES

- Baraitser, M. and Preece, M.A. 1983. The Rubinstein-Taybi syndrome: occurrence in two sets of identical twins. *Clin. Genet.* 23: 318-320.
- Shahrestanifar, M., et al. 1994. Cloning of a human cDNA encoding a putative nucleotide-binding protein related to *Escherichia coli* MinD. *Gene* 147: 281-285.
- Nakashima, H., et al. 1999. Two novel mouse genes—Nubp2, mapped to the t-complex on chromosome 17, and Nubp1, mapped to chromosome 16—establish a new gene family of nucleotide-binding proteins in eukaryotes. *Genomics* 60: 152-160.
- Bomont, P., et al. 2000. The gene encoding gigaxonin, a new member of the cytoskeletal BTB/kelch repeat family, is mutated in giant axonal neuropathy. *Nat. Genet.* 26: 370-374.
- Mathew, C.G. and Lewis, C.M. 2004. Genetics of inflammatory bowel disease: progress and prospects. *Hum. Mol. Genet.* 13: R161-R168.
- Hausmann, A., et al. 2005. The eukaryotic P loop NTPase Nbp35: an essential component of the cytosolic and nuclear iron-sulfur protein assembly machinery. *Proc. Natl. Acad. Sci. USA* 102: 3266-3271.
- Online Mendelian Inheritance in Man, OMIM™. 2007. Johns Hopkins University, Baltimore, MD. MIM Number: 600280. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Stehling, O., et al. 2008. Human Nbp35 is essential for both cytosolic iron-sulfur protein assembly and iron homeostasis. *Mol. Cell. Biol.* 28: 5517-5528.

## CHROMOSOMAL LOCATION

Genetic locus: NUBP1 (human) mapping to 16p13.13; Nubp1 (mouse) mapping to 16 A1.

## SOURCE

NUBP1 (M-13) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of NUBP1 of mouse origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-160614 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

NUBP1 (M-13) is recommended for detection of NUBP1 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with NUBP2 or NUBPL.


NUBP1 (M-13) is also recommended for detection of NUBP1 in additional species, including bovine and porcine.

Suitable for use as control antibody for NUBP1 siRNA (h): sc-93202, NUBP1 siRNA (m): sc-150090, NUBP1 shRNA Plasmid (h): sc-93202-SH, NUBP1 shRNA Plasmid (m): sc-150090-SH, NUBP1 shRNA (h) Lentiviral Particles: sc-93202-V and NUBP1 shRNA (m) Lentiviral Particles: sc-150090-V.


Molecular Weight of NUBP1: 35 kDa.

Positive Controls: COLO 320DM cell lysate: sc-2226, HeLa whole cell lysate: sc-2200 or NUBP1 (m): 293T Lysate: sc-122147.

## DATA



NUBP1 (M-13): sc-160614. Western blot analysis of NUBP1 expression in HeLa (A), COLO 320DM (B) and Jurkat (C) whole cell lysates.



NUBP1 (M-13): sc-160614. Western blot analysis of NUBP1 expression in non-transfected: sc-117752 (A) and mouse NUBP1 transfected: sc-122147 (B) 293T whole cell lysates.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

**MONOS**  
Satisfaction  
Guaranteed

Try **NUBP1 (C-3): sc-398368** or **NUBP1 (C-7): sc-514175**, our highly recommended monoclonal alternatives to NUBP1 (M-13).