

ZnT-10 (E-13): sc-160941



The Power to Question

BACKGROUND

Zinc, an essential element required for cell proliferation and differentiation, plays a role in a diverse array of cellular functions (such as neuroregulation) and acts as a cofactor for numerous enzymes and transcription factors. The zinc transporter (ZnT) family regulates the supply of zinc within cells, and its members commonly contain six membrane-spanning domains, a large histidine-rich intracellular loop and a C-terminal tail. ZnT-10 (zinc transporter 10), also known as ZNT8, ZRC1, ZNT10 or SLC30A10, is a 496 amino acid multi-pass membrane protein belonging to the cation diffusion facilitator (CDF) transporter family. Specifically expressed in fetal liver and fetal brain, ZnT-10 may be involved in the efflux of zinc from the cell. Existing as three isoforms produced by alternative splicing events, ZnT-10 enhances glucose-stimulated Insulin secretion.

REFERENCES

- Seve, M., Chimienti, F., Devergnas, S. and Favier, A. 2004. *In silico* identification and expression of SLC30 family genes: an expressed sequence tag data mining strategy for the characterization of zinc transporters' tissue expression. *BMC Genomics* 5: 32.
- Liuzzi, J.P., Bobo, J.A., Lichten, L.A., Samuelson, D.A. and Cousins, R.J. 2004. Responsive transporter genes within the murine intestinal-pancreatic axis form a basis of zinc homeostasis. *Proc. Natl. Acad. Sci. USA* 101: 14355-14360.
- Devirgiliis, C., Zalewski, P.D., Perozzi, G. and Murgia, C. 2007. Zinc fluxes and zinc transporter genes in chronic diseases. *Mutat. Res.* 622: 84-93.
- Egefjord, L., Jensen, J.L., Bang-Berthelsen, C.H., Petersen, A.B., Smidt, K., Schmitz, O., Karlsen, A.E., Pociot, F., Chimienti, F., Rungby, J. and Magnusson, N.E. 2009. Zinc transporter gene expression is regulated by pro-inflammatory cytokines: a potential role for zinc transporters in β -cell apoptosis? *BMC Endocr Disord.* 9: 7.
- Wenzlau, J.M., Frisch, L.M., Gardner, T.J., Sarkar, S., Hutton, J.C. and Davidson, H.W. 2009. Novel antigens in type 1 diabetes: the importance of ZnT8. *Curr. Diab. Rep.* 9: 105-112.
- Nicolson, T.J., Bellomo, E.A., Wijesekara, N., Loder, M.K., Baldwin, J.M., Gyulkhandanyan, A.V., Koshkin, V., Tarasov, A.I., Carzaniga, R., Kronenberger, K., Taneja, T.K., da Silva Xavier, G., Libert, S., Froguel, P., Scharfmann, R., Stetsyuk, V., Ravassard, P., Parker, H., Gribble, F.M., Reimann, F., Sladek, R., Hughes, S.J., Johnson, P.R., Masseboeuf, M., Burcelin, R., Baldwin, S.A., Liu, M., Lara-Lemus, R., Arvan, P., Schuit, F.C., Wheeler, M.B., Chimienti, F. and Rutter, G.A. 2009. Insulin storage and glucose homeostasis in mice null for the granule zinc transporter ZnT8 and studies of the type 2 diabetes-associated variants. *Diabetes* 58: 2070-2083.
- Lemaire, K., Ravier, M.A., Schraenen, A., Creemers, J.W., Van de Plas, R., Granvik, M., Van Lommel, L., Waelkens, E., Chimienti, F., Rutter, G.A., Gilon, P., in't Veld, P.A. and Schuit, F.C. 2009. Insulin crystallization depends on zinc transporter ZnT8 expression, but is not required for normal glucose homeostasis in mice. *Proc. Natl. Acad. Sci. USA* 106: 14872-14877.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

CHROMOSOMAL LOCATION

Genetic locus: SLC30A10 (human) mapping to 1q41; Slc30a10 (mouse) mapping to 1 H5.

SOURCE

ZnT-10 (E-13) is an affinity purified goat polyclonal antibody raised against a peptide mapping within a cytoplasmic domain of ZnT-10 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-160941 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

ZnT-10 (E-13) is recommended for detection of ZnT-10 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other ZnT family members.

Suitable for use as control antibody for ZnT-10 siRNA (h): sc-88233, ZnT-10 siRNA (m): sc-155818, ZnT-10 shRNA Plasmid (h): sc-88233-SH, ZnT-10 shRNA Plasmid (m): sc-155818-SH, ZnT-10 shRNA (h) Lentiviral Particles: sc-88233-V and ZnT-10 shRNA (m) Lentiviral Particles: sc-155818-V.

Molecular Weight of ZnT-10: 54 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.