

BBP (E-14): sc-161379

BACKGROUND

Proteolytic cleavage of the Amyloid protein precursor (APP) gives rise to the β -Amyloid and Amyloid A4 proteins, which are present in human platelets. Amyloid deposition is associated with type II diabetes, Down syndrome and a variety of neurological disorders, including Alzheimer's disease. Proteolytic cleavage of APP leads to the formation of the Amyloid β /A4 Amyloid protein. This protein is involved in the formation of neurofibrillary tangles and plaques that characterize the senile plaques of Alzheimer's patients. BBP (β -Amyloid-binding protein), also known as TM2D1 (TM2 domain-containing protein 1), is a 207 amino acid multi-pass membrane protein containing a G protein-coupling module that allows for interaction with the β -Amyloid peptide of APP. In cell culture, expression of BBP induces caspase-dependent vulnerability to β -Amyloid peptide toxicity, suggesting that it is a target of β -Amyloid and may be involved in the molecular pathophysiology of Alzheimer's disease.

REFERENCES

1. Kajkowski, E.M., Lo, C.F., Ning, X., Walker, S., Sofia, H.J., Wang, W., Edris, W., Chanda, P., Wagner, E., Vile, S., Ryan, K., McHendry-Rinde, B., Smith, S.C., Wood, A., Rhodes, K.J., Kennedy, J.D., Bard, J., Jacobsen, J.S. and Ozenberger, B.A. 2001. β -Amyloid peptide-induced apoptosis regulated by a novel protein containing a g protein activation module. *J. Biol. Chem.* 276: 18748-18756.
2. Taru, H., Iijima, K., Hase, M., Kirino, Y., Yagi, Y. and Suzuki, T. 2002. Interaction of Alzheimer's β -Amyloid precursor family proteins with scaffold proteins of the JNK signaling cascade. *J. Biol. Chem.* 277: 20070-20078.
3. Lee, Y., Chang, D.J., Lee, Y.S., Chang, K.A., Kim, H., Yoon, J.S., Lee, S., Suh, Y.H. and Kaang, B.K. 2003. β -Amyloid peptide binding protein does not couple to G protein in a heterologous *Xenopus* expression system. *J. Neurosci. Res.* 73: 255-259.
4. Kawasumi, M., Matsuda, S., Matsuoka, M. and Nishimoto, I. 2004. Cytoplasmic tail adaptors of Alzheimer's Amyloid- β protein precursor. *Mol. Neurobiol.* 30: 185-200.
5. Online Mendelian Inheritance in Man, OMIM™. 2006. Johns Hopkins University, Baltimore, MD. MIM Number: 610080. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
6. Cohen, L.S., Arshava, B., Estephan, R., Englander, J., Kim, H., Hauser, M., Zerbe, O., Ceruso, M., Becker, J.M. and Naider, F. 2008. Expression and biophysical analysis of two double-transmembrane domain-containing fragments from a yeast G protein-coupled receptor. *Biopolymers* 90: 117-130.
7. Chow, V.W., Mattson, M.P., Wong, P.C. and Gleichmann, M. 2010. An overview of APP processing enzymes and products. *Neuromolecular Med.* 12: 1-12.
8. Vetrivel, K.S. and Thinakaran, G. 2010. Membrane rafts in Alzheimer's disease β -Amyloid production. *Biochim. Biophys. Acta* 1801: 860-867.

CHROMOSOMAL LOCATION

Genetic locus: TM2D1 (human) mapping to 1p31.3; Tm2d1 (mouse) mapping to 4 C6.

SOURCE

BBP (E-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of BBP of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-161379 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

BBP (E-14) is recommended for detection of BBP of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

BBP (E-14) is also recommended for detection of BBP in additional species, including equine, canine and bovine.

Suitable for use as control antibody for BBP siRNA (h): sc-78572, BBP siRNA (m): sc-141482, BBP shRNA Plasmid (h): sc-78572-SH, BBP shRNA Plasmid (m): sc-141482-SH, BBP shRNA (h) Lentiviral Particles: sc-78572-V and BBP shRNA (m) Lentiviral Particles: sc-141482-V.

Molecular Weight of BBP: 19 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.