

CD1C (T-14): sc-161452



The Power to Question

BACKGROUND

The CD1 multigene family encodes five forms of the CD1 T cell surface glycoprotein in human, designated CD1A, 1B, 1C, 1D and 1E. CD1, a type 1 membrane protein, has structural similarity to the MHC class I antigen and has been shown to present lipid antigens for recognition by T lymphocytes. CD1 antigens are associated with β -2-Microglobulin and expressed on cortical thymocytes, Langerhans cells, a B cell subset and some dendritic cells.

Specifically, CD1A is a marker for Langerhans cell histiocytosis (LCH) and is found on interdigitating cells. Adaptor-protein complexes and CD1-associated chaperones control CD1 trafficking, and the development and activation of CD1-restricted T cells. Constitutive endocytosis of CD1B molecules and the differential sorting of MHC class II from lysosomes separate peptide- and lipid antigen-presenting molecules during dendritic cell maturation. CD1B is also expressed in interdigitating cells. The human CD1 genes are all closely linked in a cluster mapping at chromosome 1q 22-23.

REFERENCES

- Martin, L.H., Calabi, F., Lefebvre, F.A., Bilsland, C.A. and Milstein, C. 1987. Structure and expression of the human thymocyte antigens CD1A, CD1C, and CD1C. Proc. Natl. Acad. Sci. USA 84: 9189-9193.
- Aruffo, A. and Seed, B. 1989. Expression of cDNA clones encoding the thymocyte antigens CD1A, B, C demonstrates a hierarchy of exclusion in fibroblasts. J. Immunol. 143: 1723-1730.
- Longley, J., Kraus, J., Alonso, M. and Edelson, R. 1989. Molecular cloning of CD1A (T6), a human epidermal dendritic cell marker related to class I MHC molecules. J. Invest. Dermatol. 92: 628-631.
- Sotzki, F., Boyd, A. and Shortman, K. 1993. Surface antigens of human thymocyte populations defined by CD3, CD4 and CD8 expression: CD1A is expressed by mature thymocytes but not peripheral T cells. Immunol. Lett. 36: 101-106.
- Porcelli, S.A. 1995. The CD1 family: a third lineage of antigen-presenting molecules. Adv. Immunol. 59: 1-18.
- Melian, A., Beckman, E.M., Porcelli, S.A. and Brenner, M.B. 1996. Antigen presentation by CD1 and MHC-encoded class I-like molecules. Curr. Opin. Immunol. 8: 82-88.
- Storkus, W.J., Wei, M., Cresswell, P. and Dawson, J.R. 1996. Class I-like CD1A-C do not protect target cells from NK-mediated cytotoxicity. Cell. Immunol. 167: 154-156.
- Bauer, A., Huttinger, R., Staffler, G., Hansmann, C., Schmidt, W., Majdic, O., Knapp, W. and Stockinger, H. 1997. Analysis of the requirement for β -2-Microglobulin for expression and formation of human CD1 antigens. Eur. J. Immunol. 27: 1366-1373.
- Brigl, M. and Brenner, M.B. 2004. CD1: antigen presentation and T cell function. Annu. Rev. Immunol. 22: 817-890.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

CHROMOSOMAL LOCATION

Genetic locus: CD1C (human) mapping to 1q23.1.

SOURCE

CD1C (T-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an extracellular domain of CD1C of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-161452 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

CD1C (T-14) is recommended for detection of CD1C of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for CD1C siRNA (h): sc-42746, CD1C shRNA Plasmid (h): sc-42746-SH and CD1C shRNA (h) Lentiviral Particles: sc-42746-V.

Molecular Weight of CD1C: 43 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.