

# CDK2AP2 (G-13): sc-161474

## BACKGROUND

CDK2AP1 (cyclin-dependent kinase 2 associated protein 1), also designated p12 and DOC-1, is a growth suppressor that binds to and inhibits DNA pol  $\alpha$ /primase. When bound, CDK2AP1 affects the initiation step, but not the elongation phase, of replication. CDK2AP1 also binds to cyclin-dependent kinase 2 (Cdk2) and targets it for proteolysis. CDK2AP1 promotes cell cycle arrest by regulating the S phase of the cycle, and may trigger apoptosis. The growth factor TGF $\beta$ 1 transcriptionally-induces CDK2AP1 expression, which, in turn, mediates the growth inhibitory activity of TGF $\beta$ 1 by modulating Cdk2 activities and pRB phosphorylation. Due to its ability to trigger apoptosis, CDK2AP1 is thought to function as a tumor suppressor. CDK2AP2 (cyclin-dependent kinase 2 associated protein 2), also known as p14 or DOC-1R, is a 126 amino acid ubiquitously expressed protein that interacts with CDK2AP1 and may play a role in growth and tumor suppression.

## REFERENCES

- Zhang, X., Tsao, H., Tsuji, T., Minoshima, S., McBride, J., Majewski, P., Todd, R., Shimizu, N., Wong, D.T., Housman, D.E. and Haluska, F.G. 1999. Identification and mutation analysis of DOC-1R, a DOC-1 growth suppressor-related gene. *Biochem. Biophys. Res. Commun.* 255: 59-63.
- Tsuji, T., Duh, F.M., Latif, F., Popescu, N.C., Zimonjic, D.B., McBride, J., Matsuo, K., Ohshima, H., Todd, R., Nagata, E., Terakado, N., Sasaki, A., Matsumura, T., Lerman, M.I. and Wong, D.T. 1998. Cloning, mapping, expression, function, and mutation analyses of the human ortholog of the hamster putative tumor suppressor gene Doc-1. *J. Biol. Chem.* 273: 6704-6709.
- Shintani, S., Ohshima, H., Zhang, X., McBride, J., Matsuo, K., Tsuji, T., Hu, M.G., Hu, G., Kohno, Y., Lerman, M., Todd, R. and Wong, D.T. 2000. p12 (DOC-1) is a novel cyclin-dependent kinase 2-associated protein. *Mol. Cell Biol.* 20: 6300-6307.
- Matsuo, K., Shintani, S., Tsuji, T., Nagata, E., Lerman, M., McBride, J., Nakahara, Y., Ohshima, H., Todd, R. and Wong, D.T. 2000. p12 (DOC-1), a growth suppressor, associates with DNA polymerase  $\alpha$ /primase. *FASEB J.* 14: 1318-1324.
- Hu, M.G., Hu, G.F., Kim, Y., Tsuji, T., McBride, J., Hinds, P. and Wong, D.T. 2004. Role of p12 (CDK2AP1) in transforming growth factor- $\beta$ 1-mediated growth suppression. *Cancer Res.* 64: 490-499.
- Buajeeb, W., Zhang, X., Ohshima, H., Han, D., Surarit, R., Kim, Y. and Wong, D.T. 2004. Interaction of the Cdk2-associated protein-1, p12 (DOC-1/CDK2AP1), with its homolog, p14 (DOC-1R). *Biochem. Biophys. Res. Commun.* 315: 998-1003.
- Yuan, Z., Gaba, A.G., Kent, T.S., Bennett, A., Miller, A. and Weber, T.K. 2005. Modulation of CDK2AP1 (p12 (DOC-1)) expression in human colorectal cancer. *Oncogene* 24: 3657-3668.

## STORAGE

Store at 4°C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## CHROMOSOMAL LOCATION

Genetic locus: CDK2AP2 (human) mapping to 11q13.1; Cdk2ap2 (mouse) mapping to 19 A.

## SOURCE

CDK2AP2 (G-13) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of CDK2AP2 of human origin.

## PRODUCT

Each vial contains 200  $\mu$ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-161474 P, (100  $\mu$ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

CDK2AP2 (G-13) is recommended for detection of CDK2AP2 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with CDK2AP1.

Suitable for use as control antibody for CDK2AP2 siRNA (h): sc-96948, CDK2AP2 siRNA (m): sc-142225, CDK2AP2 shRNA Plasmid (h): sc-96948-SH, CDK2AP2 shRNA Plasmid (m): sc-142225-SH, CDK2AP2 shRNA (h) Lentiviral Particles: sc-96948-V and CDK2AP2 shRNA (m) Lentiviral Particles: sc-142225-V.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.