

MON2 (Q-13): sc-161873

BACKGROUND

MON2, also known as SF21, is a 1,718 amino acid protein that exists as multiple alternatively spliced isoforms and plays an important role in membrane trafficking. Related to the guanine nucleotide exchange factors (GEFs), MON2 shares significant homology with BIG as well as the GBF (Golgi Brefeldin A resistance factor) subfamilies of proteins. MON2 acts as a scaffold protein when associated with Dopey-1, a large cytoplasmic protein involved in trafficking between the late golgi and early endosomes. MON2 is homologous to the yeast protein and is encoded by a gene located on human chromosome 12, which encodes over 1,100 genes and comprises approximately 4.5% of the human genome. Chromosome 12 is associated with a variety of diseases and afflictions, including hypochondrogenesis, achondrogenesis, Kniest dysplasia, Noonan syndrome and trisomy 12p, which causes facial developmental defects and seizure disorders.

REFERENCES

- Allen, T.L., Brothman, A.R., Carey, J.C. and Chance, P.F. 1996. Cytogenetic and molecular analysis in trisomy 12p. *Am. J. Med. Genet.* 63: 250-256.
- Pascon, R.C. and Miller, B.L. 2000. Morphogenesis in *Aspergillus nidulans* requires Dopey (DopA), a member of a novel family of leucine zipper-like proteins conserved from yeast to humans. *Mol. Microbiol.* 36: 1250-1264.
- Avaro, S., Belgareh-Touze, N., Sibella-Argüelles, C., Volland, C. and Haguenaer-Tsapis, R. 2002. Mutants defective in secretory/vacuolar pathways in the EUROFAN collection of yeast disruptants. *Yeast* 19: 351-371.
- Chantalat, S., Courbeyrette, R., Senic-Matuglia, F., Jackson, C.L., Goud, B. and Peyroche, A. 2003. A novel Golgi membrane protein is a partner of the ARF exchange factors Gea1p and Gea2p. *Mol. Biol. Cell* 14: 2357-2371.
- Zumkeller, W., Volleth, M., Muschke, P., Tönnies, H., Heller, A., Liehr, T., Wieacker, P. and Stumm, M. 2004. Genotype/phenotype analysis in a patient with pure and complete trisomy 12p. *Am. J. Med. Genet. A* 129: 261-264.
- Efe, J.A., Plattner, F., Hulo, N., Kressler, D., Emr, S.D. and Deloche, O. 2005. Yeast Mon2p is a highly conserved protein that functions in the cytoplasm-to-vacuole transport pathway and is required for Golgi homeostasis. *J. Cell Sci.* 118: 4751-4764.
- Rachidi, M., Lopes, C., Costantine, M. and Delabar, J.M. 2005. C21orf5, a new member of Dopey family involved in morphogenesis, could participate in neurological alterations and mental retardation in Down syndrome. *DNA Res.* 12: 203-210.
- Gillingham, A.K., Whyte, J.R., Panic, B. and Munro, S. 2006. MON2, a relative of large ARF exchange factors, recruits Dop1 to the Golgi apparatus. *J. Biol. Chem.* 281 2273-2280.
- Kanjou, N., Nagao, A., Ohmiya, Y. and Ohgiya, S. 2007. Yeast mutant with efficient secretion identified by a novel secretory reporter, Cluc. *Biochem. Biophys. Res. Commun.* 358: 429-434.

CHROMOSOMAL LOCATION

Genetic locus: MON2 (human) mapping to 12q14.1; Mon2 (mouse) mapping to 10 D2.

SOURCE

MON2 (Q-13) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of MON2 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-161873 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

MON2 (Q-13) is recommended for detection of MON2 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

MON2 (Q-13) is also recommended for detection of MON2 in additional species, including equine, canine, bovine, porcine and avian.

Suitable for use as control antibody for MON2 siRNA (h): sc-96032, MON2 siRNA (m): sc-149499, MON2 shRNA Plasmid (h): sc-96032-SH, MON2 shRNA Plasmid (m): sc-149499-SH, MON2 shRNA (h) Lentiviral Particles: sc-96032-V and MON2 shRNA (m) Lentiviral Particles: sc-149499-V.

Molecular Weight of MON2: 190 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.