

TLE6 (K-14): sc-162321

BACKGROUND

The Notch signaling pathway controls various cellular interactions that are important for the specification of a variety of fates in both vertebrates and invertebrates. Key players in the Notch pathway are the TLE genes (for transducin-like enhancer of split, also designated ESG for enhancer of split groucho), which are human homologs of the *Drosophila* groucho gene. TLE6 (transducin-like enhancer of split 6), also known as GRG6, is a 449 amino acid cytoplasmic protein belonging to the WD repeat groucho/TLE family. As a member of the subcortical maternal complex (SCMC), TLE6 is essential for zygotes to progress beyond the first embryonic cell divisions. TLE6 contains seven WD repeats, a motif known to mediate protein-protein interactions. The WD40 repeat family of proteins is suggested to be involved in signal transduction, RNA processing, gene regulation, vesicular trafficking, cytoskeletal assembly and may play a role in the control of cytotypic differentiation.

REFERENCES

1. Tetsuka, T., et al. 2000. Inhibition of nuclear factor- κ B-mediated transcription by association with the amino-terminal enhancer of split, a groucho-related protein lacking WD40 repeats. *J. Biol. Chem.* 275: 4383-4390.
2. Wang, J.C., et al. 2000. Transducin-like enhancer of split proteins, the human homologs of *Drosophila* groucho, interact with hepatic nuclear factor 3 β . *J. Biol. Chem.* 275: 18418-18423.
3. Yochum, G.S. and Ayer, D.E. 2001. Pf1, a novel PHD zinc finger protein that links the TLE corepressor to the mSin3A-histone deacetylase complex. *Mol. Cell. Biol.* 21: 4110-4118.
4. Marçal, N., et al. 2005. Antagonistic effects of Grg6 and groucho/TLE on the transcription repression activity of brain factor 1/FoxG1 and cortical neuron differentiation. *Mol. Cell. Biol.* 25: 10916-10929.
5. Bajoghli, B. 2007. Evolution of the groucho/Tle gene family: gene organization and duplication events. *Dev. Genes Evol.* 217: 613-618.
6. Sekiya, T. and Zaret, K.S. 2007. Repression by groucho/TLE/Grg proteins: genomic site recruitment generates compacted chromatin *in vitro* and impairs activator binding *in vivo*. *Mol. Cell* 28: 291-303.
7. Jennings, B.H. and Ish-Horowicz, D. 2008. The groucho/TLE/Grg family of transcriptional co-repressors. *Genome Biol.* 9: 205.
8. Choi, H.K., et al. 2008. Function of multiple Lis-Homology domain/WD-40 repeat-containing proteins in feed-forward transcriptional repression by silencing mediator for retinoic and thyroid receptor/nuclear receptor corepressor complexes. *Mol. Endocrinol.* 22: 1093-1104.
9. Chen, P.C., et al. 2008. Novel roles for MLH3 deficiency and TLE6-like amplification in DNA mismatch repair-deficient gastrointestinal tumorigenesis and progression. *PLoS Genet.* 4: e1000092.

CHROMOSOMAL LOCATION

Genetic locus: TLE6 (human) mapping to 19p13.3.

SOURCE

TLE6 (K-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of TLE6 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-162321 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

TLE6 (K-14) is recommended for detection of TLE6 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other TLE family members.

TLE6 (K-14) is also recommended for detection of TLE6 in additional species, including canine.

Suitable for use as control antibody for TLE6 siRNA (h): sc-97764, TLE6 shRNA Plasmid (h): sc-97764-SH and TLE6 shRNA (h) Lentiviral Particles: sc-97764-V.

Molecular Weight of TLE6: 50 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.


 MONOS
Satisfaction
Guaranteed

Try **TLE6 (D-4): sc-515065**, our highly recommended monoclonal alternative to TLE6 (K-14).