

JRAB (E-12): sc-162967



The Power to Question

BACKGROUND

JRAB (junctional Rab13-binding protein, MICAL-like protein 2) is a 904 amino acid protein with one CH (calponin-homology) domain and one LIM zinc-binding domain. JRAB has been shown to interact with Rab13 and Rab8 to facilitate cellular transport of claudin-1, occludin and E-cadherin. This interaction is vital for the coordination of the assembly of tight junctions (TJs) and adherens junctions (AJs). Dynamic turnover (endocytic recycling) of cell-to-cell AJs and TJs is essential for epithelial morphogenesis during normal development and differentiation. The endocytic recycling of occludin and claudin proteins is part of an ongoing process of restructuring and maintaining cell junctions, especially at TJs. JRAB and Rab13 have also been implicated in the carcinoma metastasis event of epithelial cell scattering. This event shows Rab13 and JRAB colocalizing with F-actin in lamellipodial structures prior to cell scattering.

REFERENCES

1. Terai, T., Nishimura, N., Kanda, I., Yasui, N. and Sasaki, T. 2006. JRAB/MICAL-L2 is a junctional Rab13-binding protein mediating the endocytic recycling of occludin. *Mol. Biol. Cell.* 17: 2465-2475.
2. Nishimura, N. and Sasaki, T. 2008. Regulation of epithelial cell adhesion and repulsion: role of endocytic recycling. *J. Med. Invest.* 55: 9-16.
3. Nishimura, N. and Sasaki, T. 2008. Identification and characterization of JRAB/MICAL-L2, a junctional Rab13-binding protein. *Meth. Enzymol.* 438: 141-153.
4. Yamamura, R., Nishimura, N., Nakatsuji, H., Arase, S. and Sasaki, T. 2008. The interaction of JRAB/MICAL-L2 with Rab8 and Rab13 coordinates the assembly of tight junctions and adherens junctions. *Mol. Biol. Cell.* 19: 971-983.
5. Nakatsuji, H., Nishimura, N., Yamamura, R., Kanayama, H.O. and Sasaki, T. 2008. Involvement of actinin-4 in the recruitment of JRAB/MICAL-L2 to cell-cell junctions and the formation of functional tight junctions. *Mol. Cell. Biol.* 28: 3324-3335.
6. Kanda, I., Nishimura, N., Nakatsuji, H., Yamamura, R., Nakanishi, H. and Sasaki, T. 2008. Involvement of Rab13 and JRAB/MICAL-L2 in epithelial cell scattering. *Oncogene* 27: 1687-1695.
7. Nishimura, N. and Sasaki, T. 2009. Rab family small G proteins in regulation of epithelial apical junctions. *Front. Biosci.* 14: 2115-2129.
8. Sakane, A., Honda, K. and Sasaki, T. 2010. Rab13 regulates neurite outgrowth in PC12 cells through its effector protein, JRAB/MICAL-L2. *Mol. Cell. Biol.* 30: 1077-1087.

CHROMOSOMAL LOCATION

Genetic locus: MICALL2 (human) mapping to 7p22.3; Micall2 (mouse) mapping to 5 G2.

SOURCE

JRAB (E-12) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of JRAB of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-162967 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

JRAB (E-12) is recommended for detection of JRAB of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

JRAB (E-12) is also recommended for detection of JRAB in additional species, including canine and porcine.

Suitable for use as control antibody for JRAB siRNA (h): sc-89784, JRAB siRNA (m): sc-146332, JRAB shRNA Plasmid (h): sc-89784-SH, JRAB shRNA Plasmid (m): sc-146332-SH, JRAB shRNA (h) Lentiviral Particles: sc-89784-V and JRAB shRNA (m) Lentiviral Particles: sc-146332-V.

Molecular Weight of JRAB: 100 kDa.

Positive Controls: WI-38 whole cell lysate: sc-364260.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

MONOS
Satisfaction
Guaranteed

Try **JRAB (F-5): sc-376675** or **JRAB (G-10): sc-376791**, our highly recommended monoclonal alternatives to JRAB (E-12).