

α -protein kinase 3 (T-14): sc-163599

BACKGROUND

The phosphorylation and dephosphorylation of proteins on serine and threonine residues is an essential means of regulating a broad range of cellular functions in eukaryotes, including cell division, homeostasis and apoptosis. α -protein kinase 3, also known as ALPK3, MAK or KIAA1330, is a 1,907 amino acid protein that contains one α -type protein kinase domain, one Ig-like (immunoglobulin-like) domain, and belongs to the protein kinase superfamily and ALPK subfamily. Members of the ALPK subfamily include ALPK1, ALPK2 and ALPK3, all of which recognize and phosphorylate specific sites that are surrounded by peptides that have an α -helical conformation. Specifically, α -protein kinase 1 targets Myosin I and is thought to play an important role in the apical trafficking of vesicles carrying raft-associated sucrose-isomaltase (SI).

REFERENCES

- Ryazanov, A.G., Pavur, K.S. and Dorovkov, M.V. 1999. α -kinases: a new class of protein kinases with a novel catalytic domain. *Curr. Biol.* 9: R43-R45.
- Nakajima, D., Okazaki, N., Yamakawa, H., Kikuno, R., Ohara, O. and Nagase, T. 2002. Construction of expression-ready cDNA clones for KIAA genes: manual curation of 330 KIAA cDNA clones. *DNA Res.* 9: 99-106.
- Online Mendelian Inheritance in Man, OMIM[™]. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 607347. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Yamada, S., Ohira, M., Horie, H., Ando, K., Takayasu, H., Suzuki, Y., Sugano, S., Hirata, T., Goto, T., Matsunaga, T., Hiyama, E., Hayashi, Y., Ando, H., Suita, S., Kaneko, M., Sasaki, F., Hashizume, K., et al. 2004. Expression profiling and differential screening between hepatoblastomas and the corresponding normal livers: identification of high expression of the PLK1 oncogene as a poor-prognostic indicator of hepatoblastomas. *Oncogene* 23: 5901-5911.
- Drennan, D. and Ryazanov, A.G. 2004. α -kinases: analysis of the family and comparison with conventional protein kinases. *Prog. Biophys. Mol. Biol.* 85: 1-32.
- Heine, M., Cramm-Behrens, C.I., Ansari, A., Chu, H.P., Ryazanov, A.G., Naim, H.Y. and Jacob, R. 2005. α -kinase 1, a new component in apical protein transport. *J. Biol. Chem.* 280: 25637-25643.
- Thébault, S., Cao, G., Venselaar, H., Xi, Q., Bindels, R.J. and Hoenderop, J.G. 2008. Role of the α -kinase domain in transient receptor potential melastatin 6 channel and regulation by intracellular ATP. *J. Biol. Chem.* 283: 19999-20007.

CHROMOSOMAL LOCATION

Genetic locus: ALPK3 (human) mapping to 15q25.3; Alpk3 (mouse) mapping to 7 D3.

SOURCE

α -protein kinase 3 (T-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of α -protein kinase 3 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-163599 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

α -protein kinase 3 (T-14) is recommended for detection of α -protein kinase 3 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with α -protein kinase 1 or α -protein kinase 2.

α -protein kinase 3 (T-14) is also recommended for detection of α -protein kinase 3 in additional species, including equine and canine.

Suitable for use as control antibody for α -protein kinase 3 siRNA (h): sc-90044, α -protein kinase 3 siRNA (m): sc-140599, α -protein kinase 3 shRNA Plasmid (h): sc-90044-SH, α -protein kinase 3 shRNA Plasmid (m): sc-140599-SH, α -protein kinase 3 shRNA (h) Lentiviral Particles: sc-90044-V and α -protein kinase 3 shRNA (m) Lentiviral Particles: sc-140599-V.

Molecular Weight of α -protein kinase 3: 140 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker[™] compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker[™] Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz[™] Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.