

GnT-VB (Q-16): sc-164501

BACKGROUND

GnT-VB (GlcNAc-T Vb), also known as MGAT5B (mannosyl (α -1,6-)-glycoprotein β -1,6-N-acetylglucosaminyltransferase, isozyme B), GNT-IX (N-acetylglucosaminyl-transferase IX) or α -1,6-mannosylglycoprotein 6- β -N-acetylglucosaminyltransferase B, is a 792 amino acid single-pass type II membrane protein that belongs to the glycosyltransferase 18 family. Localizing to Golgi apparatus membrane, GnT-VB is predominantly expressed in brain with lower levels found in testis, spleen and thymus. GnT-VB participates in the protein glycosylation pathway and functions in the synthesis of complex cell surface N- and O-mannosyl glycans. GnT-VB also plays an important role in regulating integrin and laminin-dependent adhesion. The gene encoding GnT-VB produces five isoforms due to alternative splicing and maps to human chromosome 17q25.2 and mouse chromosome 11 E2.

REFERENCES

1. Kaneko, M., et al. 2003. A novel β (1,6)-N-acetylglucosaminyltransferase V (GnT-VB)(1). FEBS Lett. 554: 515-519.
2. Inamori, K., et al. 2003. Molecular cloning and characterization of human GnT-IX, a novel β 1,6-N-acetylglucosaminyltransferase that is specifically expressed in the brain. J. Biol. Chem. 278: 43102-43109.
3. Inamori, K., et al. 2004. N-Acetylglucosaminyltransferase IX acts on the GlcNAc β 1,2-Man α 1-Ser/Thr moiety, forming a 2,6-branched structure in brain O-mannosyl glycan. J. Biol. Chem. 279: 2337-2340.
4. Abbott, K.L., et al. 2006. Integrin-dependent neuroblastoma cell adhesion and migration on laminin is regulated by expression levels of two enzymes in the O-mannosyl-linked glycosylation pathway, PomGnT1 and GnT-Vb. Exp. Cell Res. 312: 2837-2850.
5. Lee, I., et al. 2006. N-acetylglucosaminyltransferase VB expression enhances β 1 integrin-dependent PC12 neurite outgrowth on laminin and collagen. J. Neurochem. 97: 947-956.
6. Zody, M.C., et al. 2006. DNA sequence of human chromosome 17 and analysis of rearrangement in the human lineage. Nature 440: 1045-1049.

CHROMOSOMAL LOCATION

Genetic locus: MGAT5B (human) mapping to 17q25.2; Mgat5b (mouse) mapping to 11 E2.

SOURCE

GnT-VB (Q-16) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of GnT-VB of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-164501 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

RESEARCH USE

For research use only, not for use in diagnostic procedures.

APPLICATIONS

GnT-VB (Q-16) is recommended for detection of GnT-VB of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with GnT-V.

GnT-VB (Q-16) is also recommended for detection of GnT-VB in additional species, including equine, canine and porcine.

Suitable for use as control antibody for GnT-VB siRNA (h): sc-94051, GnT-VB siRNA (m): sc-145666, GnT-VB shRNA Plasmid (h): sc-94051-SH, GnT-VB shRNA Plasmid (m): sc-145666-SH, GnT-VB shRNA (h) Lentiviral Particles: sc-94051-V and GnT-VB shRNA (m) Lentiviral Particles: sc-145666-V.


Molecular Weight of GnT-VB isoforms 1/2/3: 90/91/44 kDa.

Positive Controls: mouse brain extract: sc-2253, mouse liver extract: sc-2256 or Sol8 cell lysate: sc-2249.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



GnT-VB (Q-16): sc-164501. Western blot analysis of GnT-VB expression in mouse liver (A) and mouse brain (B) tissue extracts and Neuro-2A (C), SH-SY5Y (D) and Sol8 (E) whole cell lysates.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.