

# OTUD7A (S-15): sc-165172

## BACKGROUND

OTUD7A (OTU domain-containing protein 7A), also known as CEZANNE2, is a 926 amino acid cytoplasmic and nuclear protein that belongs to the peptidase C64 family. OTUD7A contains one A20-type zinc finger, one OTU domain and exists as two alternatively spliced isoforms. Hydrolyzing both linear and branched forms of polyubiquitin, OTUD7A has deubiquitinating activity that is directed towards "Lys-48" or "Lys-63"-linked polyubiquitin chains. The gene that encodes OTUD7A consists of approximately 174,262 bases and maps to human chromosome 15q13.3. Chromosome 15 houses over 700 genes and comprises nearly 3% of the human genome. Angelman and Prader-Willi syndromes are associated with loss of function or deletion of genes in the 15q11-q13 region. Tay-Sachs disease is a lethal disorder associated with mutations of the HEXA gene, which is encoded by chromosome 15, while Marfan syndrome is associated with chromosome 15 through the FBN1 gene.

## REFERENCES

- Knoll, J.H., Nicholls, R.D., Magenis, R.E., Graham, J.M., Lalande, M. and Latt, S.A. 1989. Angelman and Prader-Willi syndromes share a common chromosome 15 deletion but differ in parental origin of the deletion. *Am. J. Med. Genet.* 32: 285-290.
- Hurowitz, G.I., Silver, J.M., Brin, M.F., Williams, D.T. and Johnson, W.G. 1993. Neuropsychiatric aspects of adult-onset Tay-Sachs disease: two case reports with several new findings. *J. Neuropsychiatry Clin. Neurosci.* 5: 30-36.
- Boer, H., Holland, A., Whittington, J., Butler, J., Webb, T. and Clarke, D. 2002. Psychotic illness in people with Prader Willi syndrome due to chromosome 15 maternal uniparental disomy. *Lancet* 359: 135-136.
- Kayagaki, N., Phung, Q., Chan, S., Chaudhari, R., Quan, C., O'Rourke, K.M., Eby, M., Pietras, E., Cheng, G., Bazan, J.F., Zhang, Z., Arnott, D. and Dixit, V.M. 2007. DUBA: a deubiquitinase that regulates type I interferon production. *Science* 318: 1628-1632.
- Midla, G.S. 2008. Diagnosis and management of patients with Marfan syndrome. *JAAPA* 21: 21-25.
- Online Mendelian Inheritance in Man, OMIM<sup>™</sup>. 2008. Johns Hopkins University, Baltimore, MD. MIM Number: 612024. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Dan, B. 2009. Angelman syndrome: current understanding and research prospects. *Epilepsia* 50: 2331-2339.
- Ferrer-Bolufer, I., Dalmau, J., Quiroga, R., Oltra, S., Orellana, C., Monfort, S., Roselló, M., De La Osa, A. and Martínez, F. 2009. Tyrosinemia type 1 and Angelman syndrome due to paternal uniparental isodisomy 15. *J. Inherit. Metab. Dis.* 32: 349-353.

## CHROMOSOMAL LOCATION

Genetic locus: OTUD7A (human) mapping to 15q13.3; Otud7a (mouse) mapping to 7 C.

## SOURCE

OTUD7A (S-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the N-terminus of OTUD7A of human origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-165172 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

OTUD7A (S-15) is recommended for detection of OTUD7A of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with OTUD7B.

OTUD7A (S-15) is also recommended for detection of OTUD7A in additional species, including equine, canine, bovine, porcine and avian.

Suitable for use as control antibody for OTUD7A siRNA (h): sc-90013, OTUD7A siRNA (m): sc-151944, OTUD7A shRNA Plasmid (h): sc-90013-SH, OTUD7A shRNA Plasmid (m): sc-151944-SH, OTUD7A shRNA (h) Lentiviral Particles: sc-90013-V and OTUD7A shRNA (m) Lentiviral Particles: sc-151944-V.

Molecular Weight of OTUD7A: 101 kDa.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker<sup>™</sup> compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker<sup>™</sup> Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz<sup>™</sup> Mounting Medium: sc-24941.

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.