

PIG-W (P-13): sc-165248

BACKGROUND

Phosphatidylinositol-glycans (PIGs) are multi-pass transmembrane proteins that localize to the endoplasmic reticulum. PIGs exhibit various functions but all are crucial for the biosynthesis of the glycosylphosphatidylinositol (GPI)-anchor. Some PIG proteins are components of the GPI transamidase complex and play a role in the recognition of either the GPI attachment signal or the lipid portion of GPI. Other PIGs belong to the glycosyltransferase complex and function in the transfer of N-acetylglucosamine (GlcNAc) to phosphatidylinositol (PI). A variety of other PIGs play distinct roles in GPI synthesis including mannosylation of the GPI-anchor. PIG-W (phosphatidylinositol-glycan biosynthesis class W protein) is a 504 amino acid multi-pass membrane protein that functions in the third step of GPI biosynthesis and acylates the inositol ring of phosphatidylinositol.

REFERENCES

1. Stevens, V.L. and Raetz, C.R. 1991. Defective glycosyl phosphatidylinositol biosynthesis in extracts of three Thy-1 negative lymphoma cell mutants. *J. Biol. Chem.* 266: 10039-10042
2. Ryals, P.E., Pak, Y. and Thompson, G.A. 1991. Phosphatidylinositol-linked glycans and phosphatidylinositol-anchored proteins of *Tetrahymena mimbres*. *J. Biol. Chem.* 266: 15048-15053.
3. Scallan, B.J., Fung, W.J., Tsang, T.C., Li, S., Kado-Fong, H., Huang, K.S. and Kochan, J.P. 1991. Primary structure and functional activity of a phosphatidylinositol-glycan-specific phospholipase D. *Science* 252: 446-448.
4. Ueda, E., Nishimura, J., Kitani, T., Nasu, K., Kageyama, T., Kim, Y.U., Takeda, J. and Kinoshita, T. 1992. Deficient surface expression of glycosylphosphatidylinositol-anchored proteins in B cell lines established from patients with paroxysmal nocturnal hemoglobinuria. *Int. Immunol.* 4: 1263-1271.
5. Güther, M.L. and Ferguson, M.A. 1993. The microanalysis of glycosylphosphatidylinositol glycans. *Methods Mol. Biol.* 14: 99-117.
6. Watanabe, R., Inoue, N., Westfall, B., Taron, C.H., Orlean, P., Takeda, J. and Kinoshita, T. 1998. The first step of glycosylphosphatidylinositol biosynthesis is mediated by a complex of PIG-A, PIG-H, PIG-C and GPI1. *EMBO J.* 17: 877-885.
7. Murakami, Y., Siripanyapinyo, U., Hong, Y., Kang, J.Y., Ishihara, S., Nakakuma, H., Maeda, Y. and Kinoshita, T. 2003. PIG-W is critical for inositol acylation but not for flipping of glycosylphosphatidylinositol-anchor. *Mol. Biol. Cell.* 14: 4285-4295.
8. Online Mendelian Inheritance in Man, OMIM[™]. 2006. Johns Hopkins University, Baltimore, MD. MIM Number: 610275. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
9. Rodríguez, A.E., Couto, A., Echaide, I., Schnitger, L. and Florin-Christensen, M. 2010. *Babesia bovis* contains an abundant parasite-specific protein-free glycerophosphatidylinositol and the genes predicted for its assembly. *Vet. Parasitol.* 167: 227-235.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

CHROMOSOMAL LOCATION

Genetic locus: PIGW (human) mapping to 17q12; Pigw (mouse) mapping to 11 C.

SOURCE

PIG-W (P-13) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of PIG-W of mouse origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-165248 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

PIG-W (P-13) is recommended for detection of PIG-W of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other PIG family members .

Suitable for use as control antibody for PIG-W siRNA (h): sc-93698, PIG-W siRNA (m): sc-152257, PIG-W shRNA Plasmid (h): sc-93698-SH, PIG-W shRNA Plasmid (m): sc-152257-SH, PIG-W shRNA (h) Lentiviral Particles: sc-93698-V and PIG-W shRNA (m) Lentiviral Particles: sc-152257-V.

Molecular Weight of PIG-W: 57 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker[™] compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker[™] Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz[™] Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.