

Polycystin-L2 (K-15): sc-165273

BACKGROUND

Polycystin-L2, also known as PKD2L2 (polycystic kidney disease 2-like 2), is a 624 amino acid multi-pass membrane protein belonging to the Polycystin family. Polycystin-L2 contains six to eight transmembrane domains, two cation channel motifs and one channel pore calcium sodium motif. Although 48% identical and 67% similar to both Polycystin-2 and Polycystin-L, Polycystin-L2 does not possess putative EF-hand or coiled-coil domains. Existing as five alternatively spliced isoforms, Polycystin-L2 is expressed in testis, brain and kidney, with isoform 2 solely expressed in transformed lymphoblasts and isoform 3 ubiquitously expressed. Polycystin-L2 may function as a subunit of a cation channel and may play a role in fertilization. Overexpression of Polycystin-L2 increases the intracellular calcium concentration of MDCK cells, suggesting that Polycystin-L2 may be involved in the mid-late stage of spermatogenesis by way of intracellular calcium concentrations.

REFERENCES

- Veldhuisen, B., Spruit, L., Dauwerse, H.G., Breuning, M.H. and Peters, D.J. 1999. Genes homologous to the autosomal dominant polycystic kidney disease genes (PKD1 and PKD2). *Eur. J. Hum. Genet.* 7: 860-872.
- Guo, L., Schreiber, T.H., Weremowicz, S., Morton, C.C., Lee, C. and Zhou, J. 2000. Identification and characterization of a novel polycystin family member, polycystin-L2, in mouse and human: sequence, expression, alternative splicing, and chromosomal localization. *Genomics* 64: 241-251.
- Pritchard, L., Sloane-Stanley, J.A., Sharpe, J.A., Aspinwall, R., Lu, W., Buckle, V., Strmecki, L., Walker, D., Ward, C.J., Alpers, C.E., Zhou, J., Wood, W.G. and Harris, P.C. 2000. A human PKD1 transgene generates functional polycystin-1 in mice and is associated with a cystic phenotype. *Hum. Mol. Genet.* 9: 2617-2627.
- Deltas, C.C. 2001. Mutations of the human polycystic kidney disease 2 (PKD2) gene. *Hum. Mutat.* 18: 13-24.
- Stayner, C. and Zhou, J. 2001. Polycystin channels and kidney disease. *Trends Pharmacol. Sci.* 22: 543-546.
- Basora, N., Nomura, H., Berger, U.V., Stayner, C., Guo, L., Shen, X. and Zhou, J. 2002. Tissue and cellular localization of a novel polycystic kidney disease-like gene product, polycystin-L. *J. Am. Soc. Nephrol.* 13: 293-301.
- Chen, Y., Zhang, Z., Lv, X.Y., Wang, Y.D., Hu, Z.G., Sun, H., Tan, R.Z., Liu, Y.H., Bian, G.H., Xiao, Y., Li, Q.W., Yang, Q.T., Ai, J.Z., Feng, L., Yang, Y., Wei, Y.Q. and Zhou, Q. 2008. Expression of Pkd2l2 in testis is implicated in spermatogenesis. *Biol. Pharm. Bull.* 31: 1496-1500.

CHROMOSOMAL LOCATION

Genetic locus: PKD2L2 (human) mapping to 5q31.2; Pkd2l2 (mouse) mapping to 18 B1.

SOURCE

Polycystin-L2 (K-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping within a cytoplasmic domain of Polycystin-L2 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-165273 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

Polycystin-L2 (K-15) is recommended for detection of Polycystin-L2 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with Polycystin-L.

Polycystin-L2 (K-15) is also recommended for detection of Polycystin-L2 in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for Polycystin-L2 siRNA (h): sc-91892, Polycystin-L2 siRNA (m): sc-152387, Polycystin-L2 shRNA Plasmid (h): sc-91892-SH, Polycystin-L2 shRNA Plasmid (m): sc-152387-SH, Polycystin-L2 shRNA (h) Lentiviral Particles: sc-91892-V and Polycystin-L2 shRNA (m) Lentiviral Particles: sc-152387-V.

Molecular Weight of Polycystin-L2 isoforms 1/2/3/4/5: 74/70/73/72 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.