

ATP6L (S-12): sc-167165

BACKGROUND

ATP6L, also known as ATP6V0C, Vma3, ATPL, VATL or ATP6C, is a vacuolar-type H⁺-ATPase (V-ATPase). V-ATPases are multisubunit enzymes responsible for acidification of eukaryotic intracellular organelles. V-ATPases pump protons against an electrochemical gradient, while F-ATPases reverse the process, thereby synthesizing ATP. A peripheral V1 domain, which is responsible for ATP hydrolysis, and an integral V0 domain, which is responsible for proton translocation, compose V-ATPase. Nine subunits (A-H) make up the V1 domain and five subunits (a, d, c, c' and c'') make up the V0 domain. Consisting of 155 amino acids, ATP6L is a multi-pass membrane protein that makes up part of the V0 domain.

REFERENCES

- Gillespie, G.A., Somlo, S., Germino, G.G., Weinstat-Saslow, D. and Reeders, S.T. 1991. CpG island in the region of an autosomal dominant polycystic kidney disease locus defines the 5' end of a gene encoding a putative proton channel. *Proc. Natl. Acad. Sci. USA* 88: 4289-4293.
- Stevens, T.H. and Forgac, M. 1997. Structure, function and regulation of the vacuolar (H⁺)-ATPase. *Annu. Rev. Cell Dev. Biol.* 13: 779-808.
- Simckes, A.M., Swanson, S.K. and White, R.A. 2002. Chromosomal localization of three vacuolar-H⁺-ATPase 16 kDa subunit (ATP6V0C) genes in the murine genome. *Cytogenet. Genome Res.* 97: 111-115.
- Morel, N. 2003. Neurotransmitter release: the dark side of the vacuolar-H⁺-ATPase. *Biol. Cell.* 95: 453-457.
- Online Mendelian Inheritance in Man, OMIM[™]. 2003. Johns Hopkins University, Baltimore, MD. MIM Number: 108745. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Lee, I., Skinner, M.A., Guo, H.B., Sujana, A. and Pierce, M. 2004. Expression of the vacuolar H⁺-ATPase 16-kDa subunit results in the Triton X-100-insoluble aggregation of β 1 integrin and reduction of its cell surface expression. *J. Biol. Chem.* 279: 53007-53014.
- You, H., Jin, J., Shu, H., Yu, B., De Milito, A., Lozupone, F., Deng, Y., Tang, N., Yao, G., Fais, S., Gu, J. and Qin, W. 2009. Small interfering RNA targeting the subunit ATP6L of proton pump V-ATPase overcomes chemoresistance of breast cancer cells. *Cancer Lett.* 280: 110-119.

CHROMOSOMAL LOCATION

Genetic locus: ATP6V0C (human) mapping to 16p13.3; Atp6v0c (mouse) mapping to 17 A3.3.

SOURCE

ATP6L (S-12) is an affinity purified goat polyclonal antibody raised against a peptide mapping within a cytoplasmic domain of ATP6L of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-167165 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

ATP6L (S-12) is recommended for detection of ATP6L of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other ATP6 family members.

Suitable for use as control antibody for ATP6L siRNA (h): sc-93494, ATP6L siRNA (m): sc-141361, ATP6L shRNA Plasmid (h): sc-93494-SH, ATP6L shRNA Plasmid (m): sc-141361-SH, ATP6L shRNA (h) Lentiviral Particles: sc-93494-V and ATP6L shRNA (m) Lentiviral Particles: sc-141361-V.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker[™] compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker[™] Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz[™] Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.