

HLA-DP α 1 (T-12): sc-168096

BACKGROUND

Major histocompatibility complex (MHC) class II molecules destined for presentation to CD4⁺ helper T cells is determined by two key events. These events include the dissociation of class II-associated invariant chain peptides (CLIP) from an antigen binding groove in MHC class IIa/b dimers through the activity of MHC molecules HLA-DM and -DO, and subsequent peptide antigen binding. Accumulating in endosomal/lysosomal compartments and on the surface of B cells, HLA-DM and -DO molecules regulate the dissociation of CLIP and the subsequent binding of exogenous peptides to HLA class II molecules (HLA-DR, -DQ and -DP) by sustaining a conformation that favors peptide exchange. RFLP analysis of HLA-DM genes from rheumatoid arthritis (RA) patients suggests that certain polymorphisms are genetic factors for RA susceptibility. HLA-B belongs to the HLA class I heavy chain paralogs. Class I molecules play a central role in the immune system by presenting peptides derived from the endoplasmic reticulum lumen. HLA-B and C can form heterodimers consisting of a membrane anchored, heavy chain and a light chain (β -2-Microglobulin). Polymorphisms yield hundreds of HLA-B and C alleles.

REFERENCES

- Heyes, J., Austin, P., Bodmer, J., Bodmer, W., Madrigal, A., Mazzilli, M.C. and Trowsdale, J. 1986. Monoclonal antibodies to HLA-DP-transfected mouse L cells. Proc. Natl. Acad. Sci. USA 83: 3417-3421.
- Kropshofer, H., Vogt, A.B., Thery, C., Armandola, E.A., Li, B.C., Moldenhauer, G., Amigorena, S. and Hämmerling, G.J. 1998. A role for HLA-DO as a co-chaperone of HLA-DM in peptide loading of MHC class II molecules. EMBO J. 17: 2971-2981.
- Siegmund, T., Donner, H., Braun, J., Usadel, K.H. and Badenhop, K. 1999. HLA-DMA and HLA-DMB alleles in German patients with type 1 diabetes mellitus. Tissue Antigens 54: 291-294.
- Arndt, S.O., Vogt, A.B., Markovic-Plese, S., Martin, R., Moldenhauer, G., Wölpl, A., Sun, Y., Schadendorf, D., Hämmerling, G.J. and Kropshofer, H. 2000. Functional HLA-DM on the surface of B cells and immature dendritic cells. EMBO J. 19: 1241-1251.
- Brunet, A., Samaan, A., Deshaies, F., Kindt, T.J. and Thibodeau, J. 2000. Functional characterization of a lysosomal sorting motif in the cytoplasmic tail of HLA-DO. J. Biol. Chem. 275: 37062-37071.
- Doebele, R.C., Busch, R., Scott, H.M., Pashine, A. and Mellins, E.D. 2000. Determination of the HLA-DM interaction site on HLA-DR molecules. Immunity 13: 517-527.
- Louis-Pence, P., Kerlan-Candon, S., Morel, J., Combe, B., Clot, J., Pinet, V. and Eliaou, J.F. 2000. The downregulation of HLA-DM gene expression in rheumatoid arthritis is not related to their promoter polymorphism. J. Immunol. 16: 4861-4869.

CHROMOSOMAL LOCATION

Genetic locus: HLA-DPA1 (human) mapping to 6p21.32.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

SOURCE

HLA-DP α 1 (T-12) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an extracellular domain of HLA-DP α 1 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-168096 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

HLA-DP α 1 (T-12) is recommended for detection of HLA-DP α 1 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for HLA-DP α 1 siRNA (h): sc-95353, HLA-DP α 1 shRNA Plasmid (h): sc-95353-SH and HLA-DP α 1 shRNA (h) Lentiviral Particles: sc-95353-V.

Molecular Weight of HLA-DP α 1: 29 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker[™] compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker[™] Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz[™] Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.