

IGHMBP2 (S-20): sc-168144

BACKGROUND

IGHMBP2 (immunoglobulin mu binding protein 2), is also known as HCSA, HMN6, CATF1, SMARD1 or SMUBP2, is a 993 amino acid nuclear and cytoplasmic protein that is ubiquitously expressed. Belonging to the DNA2/NAM7 helicase family, IGHMBP2 is a 5' to 3' helicase that unwinds RNA and DNA duplexes in an ATP-dependent reaction. IGHMBP2 also acts as a transcriptional regulator and is necessary for transcriptional activation of the flounder liver-type antifreeze protein gene. IGHMBP2 exists as a homooligomer and is part of the cytosolic ribonucleoprotein complex. Mutations in the gene encoding IGHMBP2 are suggested to lead to distal hereditary motor neuropathy type 6 (HMN6), also known as spinal muscular atrophy distal autosomal recessive 1 (DSMA1) or spinal muscular atrophy with respiratory distress 1 (SMARD1). HMN6 is characterized by weakness and wasting of distal muscles of the anterior tibial and peroneal compartments of the legs and severe respiratory distress.

REFERENCES

1. Fukita, Y., et al. 1993. The human S μ bp-2, a DNA-binding protein specific to the single-stranded guanine-rich sequence related to the immunoglobulin μ chain switch region. *J. Biol. Chem.* 268: 17463-17470.
2. Zhang, Q., et al. 1999. Smubp-2 represses the Epstein-Barr virus lytic switch promoter. *Virology* 255: 160-170.
3. Grohmann, K., et al. 2001. Mutations in the gene encoding immunoglobulin mu-binding protein 2 cause spinal muscular atrophy with respiratory distress type 1. *Nat. Genet.* 29: 75-77.
4. Grohmann, K., et al. 2003. Infantile spinal muscular atrophy with respiratory distress type 1 (SMARD1). *Ann. Neurol.* 54: 719-724.
5. Liepinsh, E., et al. 2003. Solution structure of the R3H domain from human Smubp-2. *J. Mol. Biol.* 326: 217-223.
6. Maystadt, I., et al. 2004. Allelic heterogeneity of SMARD1 at the IGHMBP2 locus. *Hum. Mutat.* 23: 525-526.
7. Tachi, N., et al. 2005. A new mutation of IGHMBP2 gene in spinal muscular atrophy with respiratory distress type 1. *Pediatr. Neurol.* 32: 288-290.
8. de Planell-Saguer, M., et al. 2009. Biochemical and genetic evidence for a role of IGHMBP2 in the translational machinery. *Hum. Mol. Genet.* 18: 2115-2126.
9. Guenther, U.P., et al. 2009. Clinical variability in distal spinal muscular atrophy type 1 (DSMA1): determination of steady-state IGHMBP2 protein levels in five patients with infantile and juvenile disease. *J. Mol. Med.* 87: 31-41.

CHROMOSOMAL LOCATION

Genetic locus: IGHMBP2 (human) mapping to 11q13.3; Ighmbp2 (mouse) mapping to 19 A.

SOURCE

IGHMBP2 (S-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of IGHMBP2 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-168144 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

Available as TransCruz reagent for Gel Supershift and ChIP applications, sc-168144 X, 200 μ g/0.1 ml.

APPLICATIONS

IGHMBP2 (S-20) is recommended for detection of IGHMBP2 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

IGHMBP2 (S-20) is also recommended for detection of IGHMBP2 in additional species, including equine, canine and porcine.

Suitable for use as control antibody for IGHMBP2 siRNA (h): sc-96440, IGHMBP2 siRNA (m): sc-146184, IGHMBP2 shRNA Plasmid (h): sc-96440-SH, IGHMBP2 shRNA Plasmid (m): sc-146184-SH, IGHMBP2 shRNA (h) Lentiviral Particles: sc-96440-V and IGHMBP2 shRNA (m) Lentiviral Particles: sc-146184-V.

IGHMBP2 (S-20) X TransCruz antibody is recommended for Gel Supershift and ChIP applications.

Molecular Weight of IGHMBP2: 110 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.