

# p-PYK2 (Tyr 579): sc-16822

## BACKGROUND

PYK2 (proline-rich tyrosine kinase 2), a putative member of the FAK family, exhibits 61% sequence identity with FAK within its kinase domain. Like FAK, PYK2 has been shown to be a cytoplasmic protein tyrosine kinase, which is a substrate for the intrinsic protein tyrosine kinase activity of pp60Src. PYK2 (also designated CAK $\beta$  or RAFTK) is highly expressed in the central nervous system. PYK2 is rapidly phosphorylated on tyrosine residues in response to stimuli, which increases intracellular calcium levels and, in turn, activates members of the PKC family of kinases. Specifically, PYK2 is phosphorylated on Tyr 402 after stimulation with heregulin. This promotes the formation of a multiprotein complex that mediates the phosphorylation of p190 RhoGAP by Src. Activation of the PYK2 kinase leads to modulation of ion channel function and the activation of the MAPK signaling pathway. PYK2 also contains phosphorylation sites within the activation loop at Tyr 579 and Tyr 580 and on the potential GRB2-binding site at Tyr 881.

## REFERENCES

- Schaller, M.D., Borgman, C.A., Cobb, B.S., Vines, R.R., Reynolds, A.B. and Parsons, J.T. 1992. pp125 FAK, a structurally distinctive protein-tyrosine kinase associated with focal adhesions. *Proc. Natl. Acad. Sci. USA* 89: 5192-5196.
- Hanks, S.K., Calallo, M.B., Harper, M.C. and Patel, S.K. 1992. Focal adhesion protein-tyrosine kinase phosphorylated in response to cell attachment to fibronectin. *Proc. Natl. Acad. Sci. USA* 89: 8487-8491.
- Lipfert, L., Haimovich, B., Schaller, M.D., Cobb, B.S., Parsons, J.T. and Brugge, J.S. 1992. Integrin-dependent phosphorylation of the protein tyrosine kinase pp125 FAK in platelets. *J. Cell Biol.* 119: 905-912.
- Guan, J.-L., et al. 1992. Regulation of focal adhesion-associated protein tyrosine kinase by both cellular adhesion and oncogenic transformation. *Nature* 359: 690-692.

## CHROMOSOMAL LOCATION

Genetic locus: PTK2B (human) mapping to 8p21.2; Ptk2b (mouse) mapping to 14 D1.

## SOURCE

p-PYK2 (Tyr 579) is available as either goat (sc-16822) or rabbit (sc-16822-R) polyclonal affinity purified antibody raised against a short amino acid sequence containing phosphorylated Tyr 579 of PYK2 of human origin.

## PRODUCT

Each vial contains 200  $\mu$ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-16822 P, (100  $\mu$ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## APPLICATIONS

p-PYK2 (Tyr 579) is recommended for detection of Tyr 579 phosphorylated PYK2 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2  $\mu$ g per 100-500  $\mu$ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for PYK2 siRNA (h): sc-36332, PYK2 siRNA (m): sc-36333, PYK2 shRNA Plasmid (h): sc-36332-SH, PYK2 shRNA Plasmid (m): sc-36333-SH, PYK2 shRNA (h) Lentiviral Particles: sc-36332-V and PYK2 shRNA (m) Lentiviral Particles: sc-36333-V.


Molecular Weight of p-PYK2: 120 kDa.

Positive Controls: PYK2 (h): 293T Lysate: sc-115595 or Jurkat + PMA cell lysate: sc-24718.


## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: for goat primary antibody (sc-14268): use donkey anti-goat IgG-HRP: sc-2020 (range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (range: 1:2000-1:5000), for rabbit primary antibody (sc-14268-R): use goat anti-rabbit IgG-HRP: sc-2004 (range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (range: 1:2000-1:5000); Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto B Blocking Reagent: sc-2335 (use 50 mM NaF, sc-24988, as diluent), Western Blotting Luminol Reagent: sc-2048 and Lambda Phosphatase: sc-200312A. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: for goat primary antibody (sc-14268): use donkey anti-goat IgG-FITC: sc-2024 (range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (range: 1:100-1:400), for rabbit primary antibody (sc-14268-R): use goat anti-rabbit IgG-FITC: sc-2012 (range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## DATA



p-PYK2 (Tyr 579): sc-16822. Western blot analysis of PYK2 phosphorylation in non-transfected: sc-117752 (A,D) and human PYK2 transfected: sc-115595 (B) 293T whole cell lysates.



Western blot analysis of PYK2 phosphorylation in non-transfected: sc-117752 (A,D), untreated human PYK2 transfected: sc-115595 (B,E) and lambda protein phosphatase (sc-200312A) treated human PYK2 transfected: sc-115595 (C,F) 293T whole cell lysates. Antibodies tested include p-PYK2 (Tyr 579)-R: sc-16822-R (A,B,C) and PYK2 (N-19): sc-1514 (D,E,F).

## RESEARCH USE

For research use only, not for use in diagnostic procedures.