

# GABA<sub>A</sub> R $\epsilon$ (T-20): sc-16875

## BACKGROUND

GAD-65 and GAD-67, glutamate decarboxylases function to catalyze the production of GABA ( $\gamma$ -aminobutyric acid). In the central nervous system GABA functions as the main inhibitory transmitter by increasing a Cl<sup>-</sup> conductance that inhibits neuronal firing. GABA has been shown to activate both ionotropic (GABA<sub>A</sub>) and metabotropic (GABA<sub>B</sub>) receptors as well as a third class of receptors called GABA<sub>C</sub>. Both GABA<sub>A</sub> and GABA<sub>C</sub> are ligand-gated ion channels, however, they are structurally and functionally distinct. Members of the GABA<sub>A</sub> receptor family include GABA<sub>A</sub> R $\alpha$ 1-6, GABA<sub>A</sub> R  $\beta$ 1-3, GABA<sub>A</sub> R $\gamma$ 1-3, GABA<sub>A</sub> R $\delta$ , GABA<sub>A</sub> R $\epsilon$ , GABA<sub>A</sub> R $\rho$ 1 and GABA<sub>A</sub> R $\rho$ 2. The GABA<sub>B</sub> family is composed of GABA<sub>B</sub> R1 $\alpha$  and GABA<sub>B</sub> R1 $\beta$ . GABA transporters have also been identified and include GABA T-1, GABA T-2 and GABA T-3 (also designated GAT-1, -2, and -3). The GABA transporters function to terminate GABA action.

## REFERENCES

1. Nelson, H., et al. 1990. Cloning of the human brain GABA transporter. FEBS Lett. 269: 181-184.
2. Cherubini, E., et al. 1991. GABA: an excitatory transmitter in early postnatal life. Trends Neurosci. 14: 515-519.
3. Borden, L.A., et al. 1992. Molecular heterogeneity of the  $\gamma$ -aminobutyric acid (GABA) transport system. Cloning of two novel high affinity GABA transporters from rat brain. J. Biol. Chem. 267: 21098-21104.
4. Dirx, R., Jr., et al. 1995. Targeting of the 67 kDa isoform of glutamic acid decarboxylase to intracellular organelles is mediated by its interaction with the NH<sub>2</sub>-terminal region of the 65 kDa isoform of glutamic acid decarboxylase. J. Biol. Chem. 270: 2241-2246.
5. Lukasiewicz, P.D. 1996. GABA<sub>C</sub> receptors in the vertebrate retina. Mol. Neurobiol. 12: 181-194.
6. Kaupmann, K., et al. 1997. Expression cloning of GABA<sub>B</sub> receptors uncovers similarity to metabotropic glutamate receptors. Nature 386: 239-246.

## CHROMOSOMAL LOCATION

Genetic locus: GABRE (human) mapping to Xq28; Gabre (mouse) mapping to X A7.3.

## SOURCE

GABA<sub>A</sub> R $\epsilon$  (T-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of GABA<sub>A</sub> R $\epsilon$  of human origin.

## PRODUCT

Each vial contains 200  $\mu$ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-16875 P, (100  $\mu$ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## STORAGE

Store at 4°C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## APPLICATIONS

GABA<sub>A</sub> R $\epsilon$  (T-20) is recommended for detection of GABA<sub>A</sub> receptor epsilon of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for GABA<sub>A</sub> R $\epsilon$  siRNA (h): sc-42445, GABA<sub>A</sub> R $\epsilon$  siRNA (m): sc-42446, GABA<sub>A</sub> R $\epsilon$  shRNA Plasmid (h): sc-42445-SH, GABA<sub>A</sub> R $\epsilon$  shRNA Plasmid (m): sc-42446-SH, GABA<sub>A</sub> R $\epsilon$  shRNA (h) Lentiviral Particles: sc-42445-V and GABA<sub>A</sub> R $\epsilon$  shRNA (m) Lentiviral Particles: sc-42446-V.


Molecular Weight of GABA<sub>A</sub> R $\epsilon$ : 58 kDa.

Positive Controls: BC<sub>3</sub>H1 cell lysate: sc-2299, JAR cell lysate: sc-2276 or Hep G2 cell lysate: sc-2227.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## DATA



GABA<sub>A</sub> R $\epsilon$  (T-20): sc-16875. Western blot analysis of GABA<sub>A</sub> R $\epsilon$  expression in BC<sub>3</sub>H1 whole cell lysate.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

**MONOS**  
Satisfaction  
Guaranteed

Try **GABA<sub>A</sub> R $\epsilon$  (E-12): sc-271668**, our highly recommended monoclonal alternative to GABA<sub>A</sub> R $\epsilon$  (T-20).