

PRRG3 (A-12): sc-169048

BACKGROUND

PRRG3 (proline-rich γ -carboxyglutamic acid protein 3), also known as TMG3 (transmembrane γ -carboxyglutamic acid protein 3) or PRGP3, is a 231 amino acid single-pass type I membrane protein found in kidney, heart and brain. PRRG3 contains an N-terminal Gla (γ -carboxy-glutamate) domain, two PPXY motifs and is encoded by a gene that maps to human chromosome Xq28. Chromosome X consists of about 153 million base pairs and nearly 1,000 genes. Color blindness, hemophilia and Duchenne muscular dystrophy are well known X chromosome-linked conditions which affect males more frequently, as males carry a single X chromosome.

REFERENCES

1. Kulman, J.D., Harris, J.E., Xie, L. and Davie, E.W. 2001. Identification of two novel transmembrane γ -carboxyglutamic acid proteins expressed broadly in fetal and adult tissues. *Proc. Natl. Acad. Sci. USA* 98: 1370-1375.
2. Bernardino-Sgherri, J., Flagiello, D. and Dutrillaux, B. 2002. Overall DNA methylation and chromatin structure of normal and abnormal X chromosomes. *Cytogenet. Genome Res.* 99: 85-91.
3. Deeb, S.S. 2005. The molecular basis of variation in human color vision. *Clin. Genet.* 67: 369-377.
4. Online Mendelian Inheritance in Man, OMIM[™]. 2008. Johns Hopkins University, Baltimore, MD. MIM Number: 300685. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
5. Helderma-van den Enden, A.T., de Jong, R., den Dunnen, J.T., Houwing-Duistermaat, J.J., Kneppers, A.L., Ginjaar, H.B., Breuning, M.H. and Bakker, E. 2009. Recurrence risk due to germ line mosaicism: Duchenne and Becker muscular dystrophy. *Clin. Genet.* 75: 465-472.

CHROMOSOMAL LOCATION

Genetic locus: PRRG3 (human) mapping to Xq28; Prrg3 (mouse) mapping to X A7.3.

SOURCE

PRRG3 (A-12) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an extracellular domain of PRRG3 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-169048 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

APPLICATIONS

PRRG3 (A-12) is recommended for detection of PRRG3 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other PRRG family members.

PRRG3 (A-12) is also recommended for detection of PRRG3 in additional species, including equine, canine and bovine.

Suitable for use as control antibody for PRRG3 siRNA (h): sc-91294, PRRG3 siRNA (m): sc-152518, PRRG3 shRNA Plasmid (h): sc-91294-SH, PRRG3 shRNA Plasmid (m): sc-152518-SH, PRRG3 shRNA (h) Lentiviral Particles: sc-91294-V and PRRG3 shRNA (m) Lentiviral Particles: sc-152518-V.

Molecular Weight of PRRG3: 26 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker[™] compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker[™] Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz[™] Mounting Medium: sc-24941.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.