

TREML4 (I-14): sc-169672

BACKGROUND

Triggering receptor expressed on myeloid cells (TREM) family of proteins are cell surface receptors that may be involved in the innate and adaptive immune response. TREM1 may play a significant role in maintaining vascular homeostasis, coagulation, and inflammation at sites of injury through its location in platelets. Defects in the gene encoding TREM1 may be the cause of the Gray platelet syndrome, a rare inherited congenital bleeding disorder caused by a reduction or absence of the platelet α -granules in blood platelets. TREM2 is expressed in cultured B cells, as well as in T cell leukemia and monocyte leukemia, and acts as a counterreceptor for B7-H3. TREM2 interaction with B7-H3 on T-cells enhances T-cell activation. TREML4 (triggering receptor expressed on myeloid cells-like 4), also known as TL4, is a 200 amino acid secreted protein containing one Ig-like V-type (immunoglobulin-like) domain. Predominantly expressed in spleen, TREML4 is suggested to be involved in recognizing dead cells by different types of phagocytes in spleen.

REFERENCES

1. Washington, A.V., et al. 2002. Initial characterization of TREM-like transcript (TLT)-1: a putative inhibitory receptor within the TREM cluster. *Blood* 100: 3822-3824.
2. Allcock, R.J., et al. 2003. The human TREM gene cluster at 6p21.1 encodes both activating and inhibitory single IgV domain receptors and includes NKp44. *Eur. J. Immunol.* 33: 567-577.
3. Gibot, S. 2005. Clinical review: role of triggering receptor expressed on myeloid cells-1 during sepsis. *Crit. Care* 9: 485-489.
4. Lu, Y.T., et al. 2006. Preparation and characterization of monoclonal antibody against protein TREM-like transcript-1 (TLT-1). *Hybridoma* 25: 20-26.
5. King, R.G., et al. 2006. Trem-like transcript 2 is expressed on cells of the myeloid/granuloid and B lymphoid lineage and is up-regulated in response to inflammation. *J. Immunol.* 176: 6012-6021.
6. Hashiguchi, M., et al. 2008. Triggering receptor expressed on myeloid cell-like transcript 2 (TLT-2) is a counter-receptor for B7-H3 and enhances T cell responses. *Proc. Natl. Acad. Sci. USA* 105: 10495-10500.
7. Ford, J.W., et al. 2009. TREM and TREM-like receptors in inflammation and disease. *Curr. Opin. Immunol.* 21: 38-46.
8. Hemmi, H., et al. 2009. A new triggering receptor expressed on myeloid cells (Trem) family member, Trem-like 4, binds to dead cells and is a DNAX activation protein 12-linked marker for subsets of mouse macrophages and dendritic cells. *J. Immunol.* 182: 1278-1286.
9. Molloy, E.J. 2009. Triggering Receptor Expressed on Myeloid Cells (TREM) family and the application of its antagonists. *Recent Pat. Antiinfect. Drug Discov.* 4: 51-56.

CHROMOSOMAL LOCATION

Genetic locus: Trem14 (mouse) mapping to 17 C.

SOURCE

TREML4 (I-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of TREML4 of mouse origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-169672 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

TREML4 (I-14) is recommended for detection of TREML4 of mouse and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other TREML family members.

Suitable for use as control antibody for TREML4 siRNA (m): sc-154631, TREML4 shRNA Plasmid (m): sc-154631-SH and TREML4 shRNA (m) Lentiviral Particles: sc-154631-V.

Molecular Weight of TREML4: 22 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.