

BRF1 (C-17): sc-17465



BACKGROUND

Brf1, a component of the *S. Cerevisiae* RNA polymerase III transcription factor, TFIIB, functions in recruitment of RNA polymerase III to the promoter for multiple rounds of transcription. Brf1 precisely positions TFIIB on the upstream promoter-less DNA by creating stable protein-protein interactions with TATA-binding protein (TBP), another component of TFIIB. In addition, Brf1 stimulates opening of the promoter, demonstrated by the fact that a deletion in the N-terminal region of Brf1 blocks opening of the promoter past the transcriptional start site. The Ty3 retrovirus-like elements preferentially insert at the initiation sites of genes transcribed by RNA polymerase III. The Brf1 subunit mediates targeting of Ty3 elements to these genes and is necessary for their integration. Brf1 specifically binds to the C34 subunit of RNA polymerase III and the TAU135 subunit of TFIIC.

REFERENCES

1. Kassavetis, G.A., Nguyen, S.T., Kobayashi, R., Kumar, A., Geiduschek, P., Pisano, M. 1995. Cloning, expression, and function of TFC5, the gene encoding the B" component of the *Saccharomyces cerevisiae* RNA polymerase III transcription factor TFIIB. Proc. Natl. Acad. Sci. USA 92: 9786-9790.
2. Whitehall, S.K., Kassavetis, G.A., Geiduschek, E.P. 1995. The symmetry of the yeast U6 RNA gene's TATA box and the orientation of the TATA-binding protein in yeast TFIIB. Genes Dev. 9: 2974-85.
3. Ishiguro, A., Kassavetis, G.A., Geiduschek, E.P. 2002. Essential roles of Bdp1, a subunit of RNA polymerase III initiation factor TFIIB, in transcription and tRNA processing. Mol. Cell Biol. 22: 3264-75.
4. Huang, Y., McGillicuddy, E., Weindel, M., Dong, S., Maraia, R.J. 2003. The fission yeast TFIIB-related factor limits RNA polymerase III to a TATA-dependent pathway of TBP recruitment. Nucleic Acids Res. 31: 2108-16.
5. Hu P, Samudre K, Wu S, Sun Y, Hernandez N. 2004. CK2 phosphorylation of Bdp1 executes cell cycle-specific RNA polymerase III transcription repression. Mol. Cell 16: 81-92.
6. SWISS-PROT/TrEMBL (P29056). World Wide Web URL: <http://www.expasy.ch/sprot>

SOURCE

BRF1 (C-17) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of BRF1 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-17465 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

APPLICATIONS

BRF1 (C-17) is recommended for detection of BRF1 (also designated TFIIB90) of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1–2 µg per 100–500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).


Molecular Weight of BRF1: 100 kDa.

Positive Controls: K-562 nuclear extract: sc-2130 or U-937 nuclear extract: sc-2156.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/ 2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



BRF1 (C-17): sc-17465. Western blot analysis of BRF1 expression in K-562 (**A**) and U-937 (**B**) nuclear extracts.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.