SANTA CRUZ BIOTECHNOLOGY, INC.

Aldose Reductase (E-14): sc-17736



BACKGROUND

Aldose reductase (also designated AKR1B1, ALDR1, ALR2 or AR) is member of the monomeric NADPH-dependent aldoketoreductase family. Aldose reductase catalyzes the reduction of various aldehydes and has been implicated in the development of diabetic complications by catalyzing the reduction of the aldehyde form of glucose, to the corresponding sugar alcohol, sorbitol. This pathway plays a minor role in glucose metabolism in most tissues, however in diabetic hyperglycemia, cells undergoing Insulin-independent uptake of glucose accumulate significant quantities of sorbitol. The resulting hyperosmotic stress to cells may be a cause of diabetic complications such as neuropathy, retinopathy, and cataracts. Aldose reductase is very similar to human aldehyde reductase, bovine prostaglandin F synthase and to the European common frog protein, rho-crystallin.

CHROMOSOMAL LOCATION

Genetic locus: AKR1B1 (human) mapping to 7q33; Akr1b3 (mouse) mapping to 6 B1.

SOURCE

Aldose Reductase (E-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of Aldose Reductase of mouse origin.

PRODUCT

Each vial contains 200 μg lgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-17736 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

RESEARCH USE

For research use only, not for use in diagnostic procedures.

APPLICATIONS

Aldose Reductase (E-14) is recommended for detection of Aldose Reductase of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Aldose Reductase (E-14) is also recommended for detection of Aldose Reductase in additional species, including canine, bovine and porcine.

Suitable for use as control antibody for Aldose Reductase siRNA (h): sc-37119, Aldose Reductase siRNA (m): sc-29670, Aldose Reductase shRNA Plasmid (h): sc-37119-SH, Aldose Reductase shRNA Plasmid (m): sc-29670-SH, Aldose Reductase shRNA (h) Lentiviral Particles: sc-37119-V and Aldose Reductase shRNA (m) Lentiviral Particles: sc-29670-V.

Molecular Weight of Aldose Reductase: 37 kDa.

Positive Controls: JAR cell lysate: sc-2276, Caki-1 cell lysate: sc-2224 or Aldose Reductase (h6): 293T Lysate: sc-158264.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker[™] compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker[™] Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz[™] Mounting Medium: sc-24941.

DATA





Aldose Reductase (E-14): sc-17736. Western blot analysis of Aldose Reductase expression in non-transfected 2931: sc-117752 (\mathbf{A}), human Aldose Reductase transfected 2931: sc-158264 (\mathbf{B}) and Caki-1 (\mathbf{C}) whole cell lysates.

Aldose Reductase (E-14): sc-17736. Immunofluorescence staining of methanol-fixed Caki-1 cells showing cytoplasmic localization.

STORAGE

Store at 4° C, **D0 NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

MONOS Satisfation Guaranteed

Try Aldose Reductase (H-6): sc-166918 or Aldose Reductase (C-1): sc-271007, our highly recommended monoclonal alternatives to Aldose Reductase (E-14).