

β-casein (C-15): sc-17968

BACKGROUND

Milk proteins are crucial for the development of all newborn mammals and caseins constitute the major proteins in mammalian milk. β- and κ-caseins are the only caseins present in human milk. The β-casein/κ-casein ratio is higher in colostrum than in transitional and mature milk and is related to a better digestibility of colostrum casein micelles by the neonate during the first days of life. Human β-casein-encoding gene (Bca) contains a highly phosphorylated site, which is responsible for the calcium-binding capacity of β-casein. A common set of transcription factors are required for the expression of β-casein. Multiple binding sites for Stat5, C/EBPβ (CCAAT/enhancer-binding protein) and several half-sites for glucocorticoid receptor (GR) are identified in the distal human enhancer of the β-casein gene. β-casein gene transcription is regulated primarily by a composite response element (CoRE), which integrates signaling from the lactogenic hormones PRL, Insulin and hydrocortisone in mammary epithelial cells. NFκB functions as a negative regulator of β-casein gene expression during pregnancy by interfering with Stat5 tyrosine phosphorylation.

REFERENCES

- Burchell, J., Bartek, J. and Taylor-Papadimitriou, J. 1985. Production and characterization of monoclonal antibodies to human casein. A monoclonal antibody that cross-reacts with casein and α-lactalbumin. *Hybridoma* 4: 341-350.
- Lonnerdal, B., Bergstrom, S., Andersson, Y., Hjalmarsson, K., Sundqvist, A.K. and Hernell, O. 1990. Cloning and sequencing of a cDNA encoding human milk β-casein. *FEBS Lett.* 269: 153-156.
- Menon, R.S., Chang, Y.F., Jeffers, K.F., Jones, C. and Ham, R.G. 1992. Regional localization of human β-casein gene (CSN2) to 4pter-q21. *Genomics* 13: 25-26.
- Hansson, L., Edlund, A., Johansson, T., Hernell, O., Stromqvist, M., Lindquist, S., Lonnerdal, B. and Bergstrom, S. 1994. Structure of the human β-casein encoding gene. *Gene* 139: 193-199.
- Winklehner-Jennewein, P., Geymayer, S., Lechner, J., Welte, T., Hansson, L., Geley, S. and Doppler, W. 1998. A distal enhancer region in the human β-casein gene mediates the response to prolactin and glucocorticoid hormones. *Gene* 217: 127-139.
- Cuilliere, M.L., Tregoeat, V., Bene, M.C., Faure, G. and Montagne, P. 1999. Changes in the κ-casein and β-casein concentrations in human milk during lactation. *J. Clin. Lab. Anal.* 13: 213-218.
- Lykos, M.A., Fligger, J.M., Staley, M.D. and Baumrucker, C.R. 2000. Autocrine Insulin-like growth factor II inhibits β-casein mRNA expression in a mammary cell line. *J. Dairy Sci.* 83: 285-295.
- Wyszomierski, S.L. and Rosen, J.M. 2001. Cooperative effects of Stat5 (signal transducer and activator of transcription 5) and C/EBP β (CCAAT/enhancer-binding protein-β) on β-casein gene transcription are mediated by the glucocorticoid receptor. *Mol. Endocrinology* 15: 228-240.

CHROMOSOMAL LOCATION

Genetic locus: CSN2 (human) mapping to 4q21.

SOURCE

β-casein (C-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the C-terminus of β-casein of human origin.

PRODUCT

Each vial contains 200 μg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-17968 P, (100 μg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

β-casein (C-15) is recommended for detection of β-casein of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000). Suitable for use as control antibody for β-casein siRNA (h): sc-40384.

Suitable for use as control antibody for β-casein siRNA (h): sc-40384, β-casein shRNA Plasmid (h): sc-40384-SH and β-casein shRNA (h) Lentiviral Particles: sc-40384-V.

Molecular Weight of β-casein: 29 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

MONOS
Satisfaction
Guaranteed

Try **β-casein (F20.14): sc-53189**, our highly recommended monoclonal alternative to β-casein (C-15).