

# PIASx (C-20): sc-18245

## BACKGROUND

The IL-6-type family of cytokines, which includes IL-6 and a number of similar cytokines and growth factors, plays a significant role in regulating gene activation, proliferation and differentiation. Transcription factors of the STAT family are involved in IL-6 family-mediated signal transduction pathways, and upon activation undergo phosphorylation, dimerization, and translocation to the nucleus. The duration and intensity of a cell's response to cytokines can be adjusted by the effect of several regulatory mechanisms. One example involves the protein inhibitor of activated signal transducer and activator of transcription (STAT) family (PIAS family) of proteins, which act as negative regulators of STATs in cytokine signaling. PIAS proteins are able to coactivate steroid receptor-dependent transcription as well. PIASx transcript is alternatively spliced to yield two protein isoforms, PIASx- $\alpha$  and PIASx- $\beta$ , which differ in their C-terminal regions. Similar to other members of the PIAS family, the predicted PIASx proteins contain a putative zinc-binding motif and a highly acidic region.

## REFERENCES

1. Akira, S., et al. 1994. Molecular cloning of APRF, a novel IFN-stimulated gene factor 3 p91-related transcription factor involved in the gp130-mediated signaling pathway. *Cell* 77: 63-71.
2. Zhong, Z., et al. 1994. Stat3: a STAT family member activated by tyrosine phosphorylation in response to epidermal growth factor and interleukin-6. *Science* 264: 95-98.
3. Heinrich, P.C., et al. 1998. Interleukin-6-type cytokine signalling through the gp130/Jak/STAT pathway. *Biochem. J.* 334: 297-314.
4. Liu, B., et al. 1998. Inhibition of Stat1-mediated gene activation by PIAS1. *Proc. Natl. Acad. Sci. USA* 95: 10626-10631.
5. Kotaja, N., et al. 2000. ARIP3 (androgen receptor-interacting protein 3) and other PIAS (protein inhibitor of activated STAT) proteins differ in their ability to modulate steroid receptor-dependent transcriptional activation. *Mol. Endocrinol.* 14: 1986-2000.
6. Liu, B. and Shuai, K. 2001. Induction of apoptosis by protein inhibitor of activated Stat1 through c-Jun N<sub>2</sub>-terminal kinase activation. *J. Biol. Chem.* 276: 36624-36631.

## CHROMOSOMAL LOCATION

Genetic locus: PIAS2 (human) mapping to 18q21.1; Pias2 (mouse) mapping to 18 E3.

## SOURCE

PIASx (C-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of PIASx of human origin.

## PRODUCT

Each vial contains 200  $\mu$ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-18245 P, (100  $\mu$ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

PIASx (C-20) is recommended for detection of PIASx of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2  $\mu$ g per 100-500  $\mu$ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).


PIASx (C-20) is also recommended for detection of PIASx in additional species, including equine, canine, porcine and avian.

Suitable for use as control antibody for PIASx siRNA (h): sc-40849 and PIASx siRNA (m): sc-40850; and as shRNA Plasmid control antibody for PIASx shRNA Plasmid (h): sc-40849-SH and PIASx shRNA Plasmid (m): sc-40850-SH.


Molecular Weight of PIASx: 62 kDa.

Positive Controls: PIASx (m): 293 Lysate: sc-111287 or Daudi+IFN- $\alpha$  cell lysate: sc-2266.

## DATA



PIASx (C-20): sc-18245. Western blot analysis of PIASx expression in non-transfected: sc-110760 (A) and mouse PIASx transfected: sc-111287 (B) 293 whole cell lysates.



PIASx (C-20): sc-18245. Western blot analysis of PIASx expression in non-transfected: sc-117752 (A) and mouse PIASx transfected: sc-110290 (B) 293T whole cell lysates.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.

**MONOS**  
Satisfaction  
Guaranteed

Try **PIASx (D-12): sc-166494**, our highly recommended monoclonal alternative to PIASx (C-20).