


# BNP (N-10): sc-18813



## BACKGROUND

Natriuretic peptides comprise a family of three structurally related molecules: atrial natriuretic peptide (ANP), brain natriuretic peptide (BNP) and C-type natriuretic peptide (CNP). ANP and BNP act mainly as cardiac hormones, produced primarily by the atrium and ventricle, respectively, while the gene encoding C-type natriuretic peptide is expressed mainly in the brain. These peptides possess potent natriuretic, diuretic and vasodilating activities and are implicated in body fluid homeostasis and blood pressure control. ANP, BNP and CNP are highly homologous within the 17-residue ring structure formed by an intramolecular disulfide linkage. The genes which encode for ANP and BNP map to human chromosome 1p36.22. The gene which encodes for CNP maps to human chromosome 2q24-qter.

## REFERENCES

1. Saito, T., Matsuda, M., Yamaguchi, T., Akatsuka, S. and Nakagawa, S. 1975. Proceedings: Systemic-pulmonary arteriovenous fistula—a report of a case. *Jpn. Circ. J.* 39: 723.
2. Mair, J., Hammerer-Lercher, A. and Puschendorf, B. 2001. The impact of cardiac natriuretic peptide determination on the diagnosis and management of heart failure. *Clin. Chem. Lab. Med.* 39: 571-588.
3. Cowie, M.R. and Mendez, G.F. 2002. BNP and congestive heart failure. *Prog. Cardiovasc. Dis.* 44: 293-321.
4. Hobbs, F.D., Davis, R.C., Roalfe, A.K., Hare, R., Davies, M.K. and Kenkre, J.E. 2002. Reliability of N-terminal pro-brain natriuretic peptide assay in diagnosis of heart failure: cohort study in representative and high risk community populations. *BMJ* 324: 1498.
5. Hall, C. 2004. Essential biochemistry and physiology of (NT-pro)BNP. *Eur. J. Heart Fail.* 6: 257-260.
6. Pfister, R., Scholz, M., Wielckens, K., Erdmann, E. and Schneider, C.A. 2004. Use of NT-proBNP in routine testing and comparison to BNP. *Eur. J. Heart Fail.* 6: 289-293.

## CHROMOSOMAL LOCATION

Genetic locus: NPPB (human) mapping to 1p36.2.

## SOURCE

BNP (N-10) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of BNP of human origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-18813 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## APPLICATIONS

BNP (N-10) is recommended for detection of precursor and mature γ-BNP, BNP 32, and natriuretic peptide of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for BNP siRNA (h): sc-43636, BNP shRNA Plasmid (h): sc-43636-SH and BNP shRNA (h) Lentiviral Particles: sc-43636-V.

Molecular Weight of glycosylated BNP precursor: 25-36 kDa.

Molecular Weight of deglycosylated mature BNP: 12 kDa.

Positive Controls: SK-N-MC cell lysate: sc-2237, IMR-32 cell lysate: sc-2409 or U87-MG cell lysate.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.



Try **BNP (D-8): sc-271185** or **BNP (8D5B4C11): sc-101418**, our highly recommended monoclonal alternatives to BNP (N-10).