

BAFF-R (M-20): sc-19889

BACKGROUND

TNF cytokine family member TALL-1 (also designated BAFF, zTNF4, THANK and BLYS) is a type II membrane protein that shares characteristics with other members of the tumor necrosis factor (TNF) cytokine family. TALL-1 has the ability to bind to three receptors, TACI, BCMA and BAFF-R, but unlike other TNF receptors, BAFF-R specifically binds only the TALL-1 ligand. The gene encoding human BAFF-R, which maps to chromosome 22q13.1-q13.31, is expressed at high levels in spleen and lymph nodes and at lower levels in peripheral blood leukocytes and thymus. Expression of BAFF-R is crucial for selecting transitional B cells into the mature B cell pool.

REFERENCES

1. Schneider, P., MacKay, F., Steiner, V., Hofmann, K., Bodmer, J.L., Holler, N., Ambrose, C., Lawton, P., Bixler, S., Acha-Orbea, H., Valmori, D., Romero, P., Werner-Favre, C., Zubler, R.H., Browning, J.L. and Tschoop, J. 1999. BAFF, a novel ligand of the tumor necrosis factor family, stimulates B cell growth. *J. Exp. Med.* 189: 1747-1756.
2. Moore, P.A., Belvedere, O., Orr, A., Pieri, K., LaFleur, D.W., Feng, P., Soppet, D., Charters, M., Gentz, R., Parmelee, D., Li, Y., Galperina, O., Giri, J., Roschke, V., Nardelli, B., Carrell, J., Sosnovtseva, S., Greenfield, W., Ruben, S.M., Olsen, H.S., Fikes, J. and Hilbert, D.M. 1999. BLYS: member of the tumor necrosis factor family and B lymphocyte stimulator. *Science* 285: 260-263.
3. Mukhopadhyay, A., Ni, J., Zhai, Y., Yu, G.L. and Aggarwal, B.B. 1999. Identification and characterization of a novel cytokine, THANK, a TNF homologue that activates apoptosis, nuclear factor-B, and c-Jun NH₂-terminal kinase. *J. Biol. Chem.* 274: 15978-15981.
4. Waldschmidt, T.J. and Noelle, R.J. 2001. Immunology. Long live the mature B cell—a baffling mystery resolved. *Science* 293: 2012-2013.
5. Thompson, J.S., Bixler, S.A., Qian, F., Vora, K., Scott, M.L., Cachero, T.G., Hession, C., Schneider, P., Sizing, I.D., Mullen, C., Strauch, K., Zafari, M., Benjamin, C.D., Tschoop, J., Browning, J.L. and Ambrose, C. 2001. BAFF-R, a newly identified TNF receptor that specifically interacts with BAFF. *Science* 293: 2108-2111.

CHROMOSOMAL LOCATION

Genetic locus: Tnfrsf13c (mouse) mapping to 15 E1.

SOURCE

BAFF-R (M-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of BAFF-R of mouse origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-19889 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

BAFF-R (M-20) is recommended for detection of BAFF-R of mouse and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for BAFF-R siRNA (m): sc-40232, BAFF-R shRNA Plasmid (m): sc-40232-SH and BAFF-R shRNA (m) Lentiviral Particles: sc-40232-V.

Molecular Weight of BAFF-R monomer: 19 kDa.

Molecular Weight of BAFF-R dimer: 40 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.