



ND5 (C-13): sc-20505

BACKGROUND

NADH-ubiquinone oxidoreductase (complex I) is a complicated multiprotein complex located in the inner mitochondrial membrane. Human complex I is important for energy metabolism because its main function is to transport electrons from NADH to ubiquinone, which is accompanied by translocation of protons from the mitochondrial matrix to the intermembrane space. Human complex I appears to consist of 41 subunits. A small number of complex I subunits are the products of mitochondrial genes (subunits 1–7), while the remainder are nuclear encoded and imported from the cytoplasm. Various tissues from patients with neurodegenerative disease are frequently deficient in complex I. The transcript expression of the complex I subunit 4 is significantly decreased in the cell models of Parkinson's disease (PD), suggesting that functional defects in complex I enzyme activity in PD may result from changes in mRNA levels. Complex I 24 kDa subunit is significantly reduced in occipital cortex and thalamus in patients with Down syndrome (DS) and temporal and occipital cortices in patients with Alzheimer's diseases (AD). Complex I 75 kDa subunit is significantly reduced in brain regions from patients with DS (temporal, occipital and caudate nucleus) and AD (parietal cortex). Thus, reductions of two subunits of complex I may lead to the impairment of energy metabolism and may result in neuronal cell death.

REFERENCES

1. Ton, C., Hwang, D.M., Dempsey, A.A., and Liew, C.C. 1997. Identification and primary structure of five human NADH-ubiquinone oxidoreductase subunits. *Biochem. Biophys. Res. Commun.* 241: 589-594.
2. Loeffen, J.L., Triepels, R.H., van den Heuvel, L.P., Schuelke, M., Buskens, C.A., Smeets, R.J., Trijbels, J.M., and Smeitink, J.A. 1998. cDNA of eight nuclear encoded subunits of NADH:ubiquinone oxidoreductase: human complex I cDNA characterization completed. *Biochem. Biophys. Res. Commun.* 253: 415-422.
3. Smeitink, J., Loeffen, J., Smeets, R., Triepels, R., Ruitenbeek, W., Trijbels, F., and van den Heuvel, L. 1998. Molecular characterization and mutational analysis of the human B17 subunit of the mitochondrial respiratory chain complex I. *Hum. Genet.* 103: 245-250.
4. Conn, K.J., Ullman, M.D., Eisenhauer, P.B., Fine, R.E., and Wells, J.M. 2001. Decreased expression of the NADH:ubiquinone oxidoreductase (complex I) subunit 4 in 1-methyl-4-phenylpyridinium-treated human neuroblastoma SH-SY5Y cells. *Neurosci. Letts.* 306: 145-148.
5. Kim, S.H., Vlkolinsky, R., Cairns, N., Fountoulakis, M., and Lubec, G. 2001. The reduction of NADH ubiquinone oxidoreductase 24- and 75-kDa subunits in brains of patients with Down syndrome and Alzheimer's disease. *Life Sci.* 68: 2741-2750.

SOURCE

ND5 (C-13) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of ND5 of human origin.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-20505 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

ND5 (C-13) is recommended for detection of NADH-ubiquinone oxidoreductase chain 5 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Molecular Weight of ND5: 70 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.