

Trk C (I-20): sc-20546

BACKGROUND

The Trk proto-oncogene encodes a tyrosine protein kinase, Trk A, also designated Trk gp140, that serves as a receptor for certain neurotrophic factors including nerve growth factor (NGF) and neurotrophin-3 (NT-3). Trk B is a tyrosine kinase gene highly related to Trk A. Trk B expression is confined to tissues within the central and peripheral nervous systems. The brain-derived neurotrophic factor (BDNF) and NT-3, but not NGF, can induce rapid phosphorylation on Tyrosine of Trk B gp145, one of the receptors encoded by Trk B, although BDNF elicits a response at least two orders of magnitude greater than NT-3. Thus it appears that Trk B gp145 may represent a neurotrophic receptor for BDNF and NT-3. The third member of the Trk family of tyrosine kinases, Trk C, is mainly expressed in nervous tissue, and is a functional receptor for NT-3. Four forms of Trk C are produced by alternative splicing. Isoform A is full length Trk C, isoform B differs from A between amino acids 529 to 612 and is missing amino acids 613 to 839. Trk C, isoform C is missing amino acids 712 to 725 and isoform D is missing amino acids 402 to 410.

REFERENCES

1. Klein, R., Parada, L.F., Coulier, F. and Barbacid, M. 1989. Trk B, a novel tyrosine protein kinase receptor expressed during mouse neural development. *EMBO J.* 8: 3701-3709.
2. Klein, R., Martin-Zanca, D., Barbacid, M. and Parada, L.F. 1990. Expression of the tyrosine kinase receptor gene Trk B is confined to the murine embryonic and adult nervous system. *Development* 109: 845-850.
3. Kaplan, D.R., Martin-Zanca, D. and Parada, L.F. 1991. Tyrosine phosphorylation and tyrosine kinase activity of the Trk proto-oncogene product induced by NGF. *Nature* 350: 158-160.
4. Hempstead, B.L., Martin-Zanca, D., Kaplan, D.R., Parada, L.F. and Chao, M.V. 1991. High-affinity NGF binding requires coexpression of the Trk proto-oncogene and the low-affinity NGF receptor. *Nature* 350: 678-683.
5. Klein, R., Jing, S., Nanduri, V., O'Rourke, E. and Barbacid, M. 1991. The Trk proto-oncogene encodes a receptor for nerve growth factor. *Cell* 65: 189-197.
6. Cordon-Cardo, C., Tapley, P., Jing, S., Nanduri, V., O'Rourke, E., Lamballe, F., Kovary, K., Klein, R., Jones, K.R., Reichardt, L.F. and Barbacid, M. 1991. The Trk tyrosine protein kinase mediates the mitogenic properties of nerve growth factor and neurotrophin-3. *Cell* 66: 173-183.

CHROMOSOMAL LOCATION

Genetic locus: NTRK3 (human) mapping to 15q25.3; Ntrk3 (mouse) mapping to 7 D3.

SOURCE

Trk C (I-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the N-terminus of Trk C of human origin.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-20546 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

Trk C (I-20) is recommended for detection of Trk C of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Trk C (I-20) is also recommended for detection of Trk C in additional species, including bovine and porcine.

Suitable for use as control antibody for Trk C siRNA (h): sc-36730, Trk C siRNA (m): sc-36731, Trk C shRNA Plasmid (h): sc-36730-SH, Trk C shRNA Plasmid (m): sc-36731-SH, Trk C shRNA (h) Lentiviral Particles: sc-36730-V and Trk C shRNA (m) Lentiviral Particles: sc-36731-V.

Molecular Weight of Trk C: 145 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.


 MONOS
 Satisfaction
 Guaranteed

Try **Trk C (WW6): sc-80403**, our highly recommended monoclonal alternative to Trk C (I-20).