

PIGF (H-90): sc-20714

BACKGROUND

The onset of angiogenesis is believed to be an early event in tumorigenesis and may facilitate tumor progression and metastasis. Several growth factors with angiogenic activity have been described. These include fibroblast growth factor (FGF), platelet derived growth factor (PDGF), vascular endothelial growth factor (VEGF) and placenta growth factor (PIGF). Like VEGF, several PIGF variants have been shown to arise from alternative mRNA splicings. Evidence has suggested VEGF to be an obligatory component in PIGF signaling. While VEGF homodimers and VEGF/PIGF heterodimers function as potent mediators of mitogenic and chemotactic responses in endothelial cells, PIGF homodimers are effectual only at extremely high concentrations. Indeed, many of the physiological effects attributed to VEGF may actually be a result of VEGF/PIGF. VEGF and PIGF share a common receptor, Flt-1, and may also activate Flk-1/KDR.

REFERENCES

1. Folkman, J., et al. 1987. Angiogenic factors. *Science* 235: 442-447.
2. Folkman, J., et al. 1989. Induction of angiogenesis during the transition from hyperplasia to neoplasia. *Nature* 339: 58-61.
3. Bouck, N. 1990. Tumor angiogenesis: the role of oncogenes and tumor suppressor genes. *Cancer Cells* 2: 179-185.
4. Ferrara, N., et al. 1991. The vascular endothelial growth factor family of polypeptides. *J. Cell. Biochem.* 47: 211-218.
5. DiSalvo, J., et al. 1995. Purification and characterization of a naturally occurring vascular endothelial growth factor placenta growth factor heterodimer. *J. Biol. Chem.* 270: 7717-7723.
6. Cao, Y., et al. 1996. Heterodimers of placenta growth factor/vascular endothelial growth factor. Endothelial activity, tumor cell expression and high affinity binding to Flk-1/KDR. *J. Biol. Chem.* 271: 3154-3162.
7. Clauss, M., et al. 1996. The vascular endothelial growth factor receptor Flt-1 mediates biological activities. Implications for functional role of placenta growth factor in monocyte activation and chemotaxis. *J. Biol. Chem.* 271: 17629-17634.

CHROMOSOMAL LOCATION

Genetic locus: PGF (human) mapping to 14q24.3.

SOURCE

PIGF (H-90) is a rabbit polyclonal antibody raised against amino acids 51-100 mapping within an internal region of PIGF isoform 1 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

APPLICATIONS

PIGF (H-90) is recommended for detection of PIGF-1, PIGF-2, and PIGF-3 isoforms of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

PIGF (H-90) is also recommended for detection of PIGF-1, PIGF-2, and PIGF-3 isoforms in additional species, including equine, canine and porcine.

Suitable for use as control antibody for PIGF siRNA (h): sc-44027, PIGF shRNA Plasmid (h): sc-44027-SH and PIGF shRNA (h) Lentiviral Particles: sc-44027-V.

Molecular Weight of PIGF: 18 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.


 MONOS
 Satisfation
 Guaranteed

Try **PIGF (MM0010-2D93): sc-101572**, our highly recommended monoclonal alternative to PIGF (H-90).