

NCCT (C-16): sc-21554



The Power to Question

BACKGROUND

NCCT (Na-Cl cotransporter), also known as NCC, TSC (thiazide-sensitive sodium-chloride cotransporter) or SLC12A3 (solute carrier family 12 (sodium/chloride transporters), member 3), is a 1,021 amino acid cell membrane protein that is predominately expressed in kidney. Belonging to the SLC12A transporter family, NCCT plays an important role in renal sodium reabsorption, functioning as a renal thiazide-sensitive sodium-chloride cotransporter. Mutations to the NCCT gene are the cause of Gitelman syndrome, an autosomal recessive disorder characterized by hypokalemic alkalosis combined with hypomagnesemia, low urinary calcium, and increased renin activity associated with normal blood pressure.

REFERENCES

1. Mastroianni, N., et al. 1996. Novel molecular variants of the Na-Cl cotransporter gene are responsible for Gitelman syndrome. *Am. J. Hum. Genet.* 59: 1019-1026.
2. Mastroianni, N., et al. 1996. Molecular cloning, expression pattern, and chromosomal localization of the human Na-Cl thiazide-sensitive cotransporter (SLC12A3). *Genomics* 35: 486-493.
3. Takeuchi, K., et al. 1996. Association of a mutation in thiazide-sensitive Na-Cl cotransporter with familial Gitelman's syndrome. *J. Clin. Endocrinol. Metab.* 81: 4496-4499.
4. Simon, D.B., et al. 1996. Gitelman's variant of Bartter's syndrome, inherited hypokalaemic alkalosis, is caused by mutations in the thiazide-sensitive Na-Cl cotransporter. *Nat. Genet.* 12: 24-30.
5. Lemmink, H.H., et al. 1998. Novel mutations in the thiazide-sensitive NaCl cotransporter gene in patients with Gitelman syndrome with predominant localization to the C-terminal domain. *Kidney Int.* 54: 720-730.
6. Melander, O., et al. 2000. Genetic variants of thiazide-sensitive NaCl cotransporter in Gitelman's syndrome and primary hypertension. *Hypertension* 36: 389-394.
7. Monkawa, T., et al. 2000. Novel mutations in thiazide-sensitive Na-Cl cotransporter gene of patients with Gitelman's syndrome. *J. Am. Soc. Nephrol.* 11: 65-70.

CHROMOSOMAL LOCATION

Genetic locus: SLC12A3 (human) mapping to 16q13; Slc12a3 (mouse) mapping to 8 C5.

SOURCE

NCCT (C-16) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the C-terminus of NCCT of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-21554 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

NCCT (C-16) is recommended for detection of NCCT of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

NCCT (C-16) is also recommended for detection of NCCT in additional species, including equine, canine, bovine and avian.

Suitable for use as control antibody for NCCT siRNA (h): sc-42515, NCCT siRNA (m): sc-42516, NCCT shRNA Plasmid (h): sc-42515-SH, NCCT shRNA Plasmid (m): sc-42516-SH, NCCT shRNA (h) Lentiviral Particles: sc-42515-V and NCCT shRNA (m) Lentiviral Particles: sc-42516-V.


Molecular Weight of NCCT: 114 kDa.

Positive Controls: KNRK whole cell lysate: sc-2214.


RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



NCCT (C-16): sc-21554. Western blot analysis of NCCT expression in KNRK whole cell lysate.



NCCT (C-16): sc-21554. Immunofluorescence staining of formalin-fixed, paraffin-embedded mouse kidney tissue showing apical membrane staining of cells in distal tubules and nuclear DAPI counterstain, at high magnification. Kindly provided by Vijay Marati Radhakrishnan, Ph.D., Steele Children's Research Center, The University of Arizona Department of Pediatrics.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.