SANTA CRUZ BIOTECHNOLOGY, INC.

ABCG8 (N-20): sc-21842



BACKGROUND

The ABCG (White) subfamily of ABC transporters, which includes ABCG1, ABCG5 and ABCG8, are critically involved in the regulation of lipid-trafficking mechanisms in macrophages, hepatocyes and intestinal mucosa cells. ABCG8 (also designated Sterolin-2) is expressed in the liver, small intestine and colon. ABCG8 normally cooperates with ABCG5 to limit intestinal absorption and to promote biliary excretion of sterols, whereas mutated forms of ABCG8 and ABCG5 cause sterol accumulation and atherosclerosis. ABCG8 and ABCG5 genes are also distinct targets of the LXR α and LXR β oxysterol receptors, which serve as sterol sensors to coordinately regulate sterol catabolism, storage, efflux and elimination. Mutations in either ABCG8 or ABCG5 lead to sitosterolemia, a rare autosomal recessive disorder characterized by hyperabsorption of all sterols, including cholesterol and plant and shellfish sterols. Patients with this disease are hypercholesterolemic and frequently develop xanthomas, accelerated atherosclerosis and premature coronary artery disease.

REFERENCES

- Berge, K.E., Tian, H., Graf, G.A., Yu, L., Grishin, N.V., Schultz, J., Kwiterovich, P., Shan, B., Barnes, R. and Hobbs, H.H. 2000. Accumulation of dietary cholesterol in sitosterolemia caused by mutations in adjacent ABC transporters. Science 290: 1771-1775.
- 2. Lee, M.H., Lu, K. and Patel, S.B. 2001. Genetic basis of sitosterolemia. Curr. Opin. Lipidol. 12: 141-149.
- Lee, M.H., Lu, K., Hazard, S., Yu, H., Shulenin, S., Hidaka, H., Kojima, H., Allikmets, R., Sakuma, N., Pegoraro, R., Srivastava, A.K., Salen, G., Dean, M. and Patel, S.B. 2001. Identification of a gene, ABCG5, important in the regulation of dietary cholesterol absorption. Nat. Genet. 27: 79-83.
- Lu, K., Lee, M.H., Hazard, S., Brooks-Wilson, A., Hidaka, H., Kojima, H., Ose, L., Stalenhoef, A.F., Mietinnen, T., Bjorkhem, I., Bruckert, E., Pandya, A., Brewer H.B., Jr., Salen, G., Dean, M., Srivastava, A. and Patel, S.B. 2001. Two genes that map to the STSL locus cause sitosterolemia: genomic structure and spectrum of mutations involving Sterolin-1 and Sterolin-2, encoded by ABCG5 and ABCG8, respectively. Am. J. Hum. Genet. 69: 278-290.
- Schmitz, G., Langmann, T. and Heimerl, S. 2001. Role of ABCG1 and other ABCG family members in lipid metabolism. J. Lipid Res. 42: 1513-1520.
- 6. Repa, J.J., Berge, K.E., Pomajzl, C., Richardson, J.A., Hobbs, H. and Mangelsdorf, D.J. 2002. Regulation of ATP-binding cassette sterol transporters ABCG5 and ABCG8 by the liver X receptors α and β . J. Biol. Chem. 277: 18793-18800.
- Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 605460. World Wide Web URL: http://www.ncbi.nlm.nih. gov/omim/

CHROMOSOMAL LOCATION

Genetic locus: ABCG8 (human) mapping to 2p21.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

SOURCE

ABCG8 (N-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the N-terminus of ABCG8 of human origin.

PRODUCT

Each vial contains 200 μg lgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-21842 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

ABCG8 (N-20) is recommended for detection of ABCG8 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for ABCG8 siRNA (h): sc-41154, ABCG8 shRNA Plasmid (h): sc-41154-SH and ABCG8 shRNA (h) Lentiviral Particles: sc-41154-V.

Molecular Weight of ABCG8: 75 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluo-rescence: use donkey anti-goat IgG-FITC: sc-2783 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **D0 NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.