

# p-Raf-1 (Ser 621)-R: sc-21874-R

## BACKGROUND

Raf-1 is an ubiquitously expressed cytoplasmic protein with intrinsic serine/threonine kinase activity. Raf-1, or c-Raf, is the cellular homolog of v-Raf, the product of the transforming gene of the 3611 strain of murine sarcoma virus. The unregulated kinase activity of the v-Raf protein is associated with cellular transformation and mitogenesis. Raf-1 is normally suppressed by its regulatory N-terminal domain. Raf-1 is activated in response to a variety of tyrosine kinase receptors as well as in response to pp60v-Src expression. Specifically, Raf-1 is phosphorylated in the catalytic domain at Ser 338 and, to a lesser extent, Ser 339. This phosphorylation requires the co-activation of PI 3-kinase and the Ras signaling pathway. Raf-1 is also phosphorylated on Tyr 340 and 341, which induces the phosphorylation of MEK. Phosphorylation of Ser 621 is essential for the catalytic activity of Raf-1 and downregulation by c-AMP-dependent protein kinase A (PKA). PKA also phosphorylates Raf-1 on Ser 43 and Ser 259. PKA phosphorylation of Ser 259 inhibits Raf-1 and decreases the phosphorylation necessary for Raf-1 activation at Ser 338.

## REFERENCES

- Rapp, U.R., Goldsborough, M.D., Mark, G.E., Bonner, T.I., Groffen, J., Reynolds, F.H. and Stephenson, J.R. 1983. Structure and biological activation of v-Raf, a unique oncogene transduced by a retrovirus. *Proc. Natl. Acad. Sci. USA* 80: 4218-4222.
- Huleihel, M., Goldsborough, M., Cleveland, J., Gunnell, M., Bonner, T. and Rapp, U.R. 1986. Characterization of murine A-Raf, a new oncogene related to the v-Raf oncogene. *Mol. Cell. Biol.* 6: 2655-2662.
- Heidecker, G., Huleihel, M., Cleveland, J.L., Kolch, W., Beck, P., Lloyd, P., Pawson, T. and Rapp, U.R. 1990. Mutational activation of c-Raf-1 and definition of the minimal transforming sequence. *Mol. Cell. Biol.* 10: 2503-2512.
- Mischak, H., Seitz, T., Janosch, P., Eulitz, M., Steen, H., Schellerer, M., Philipp, A. and Kolch W. 1996. Negative regulation of Raf-1 by phosphorylation of Serine 621. *Mol. Cell. Biol.* 16: 5409-5418.
- Diaz, B., Barnard, D., Filson, A., MacDonald, S., King, A. and Marshall, M. 1997. Phosphorylation of Raf-1 Serine 338/Serine 339 is an essential regulatory event for Ras-dependent activation and biological signaling. *Mol. Cell. Biol.* 17: 4509-4516.
- King, A.J., Sun, H., Diaz, B., Barnard, D., Miao, W., Bagrodia, S. and Marshall, M.S. 1998. The protein kinase PAK3 positively regulates Raf-1 activity through phosphorylation of Serine 338. *Nature* 396: 180-183.
- Barnard, D., Diaz, B., Clawson, D. and Marshall, M. 1998. Oncogenes, growth factors and phorbol esters regulate Raf-1 through common mechanisms. *Oncogene* 17: 1539-1547.
- Chaudhary, A., King, W.G., Mattaliano, M.D., Frost, J.A., Diaz, B., Morrison, D.K., Cobb, M.H., Marshall, M.S. and Brugge, J.S. 2000. Phosphatidylinositol 3-kinase regulates Raf-1 through PAK phosphorylation of Serine 338. *Curr. Biol.* 10: 551-554.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## CHROMOSOMAL LOCATION

Genetic locus: RAF1 (human) mapping to 3p25.2; Raf1 (mouse) mapping to 6 E3.

## SOURCE

p-Raf-1 (Ser 621)-R is a rabbit polyclonal antibody raised against a short amino acid sequence containing Ser 621 phosphorylated Raf-1 of human origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-21874 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

p-Raf-1 (Ser 621)-R is recommended for detection of Ser 621 phosphorylated Raf-1 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

p-Raf-1 (Ser 621)-R is also recommended for detection of correspondingly phosphorylated Raf-1 in additional species, including equine, canine, bovine, porcine and avian.

Suitable for use as control antibody for Raf-1 siRNA (h): sc-29462, Raf-1 siRNA (m): sc-29463, Raf-1 shRNA Plasmid (h): sc-29462-SH, Raf-1 shRNA Plasmid (m): sc-29463-SH, Raf-1 shRNA (h) Lentiviral Particles: sc-29462-V and Raf-1 shRNA (m) Lentiviral Particles: sc-29463-V.

Molecular Weight of p-Raf-1: 74 kDa.

Positive Controls: HeLa+UV irradiated cell lysate: sc-2221, HeLa + UV irradiated cell lysate: sc-2221 or HeLa whole cell lysate: sc-2200.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto B Blocking Reagent: sc-2335 (use 50 mM NaF, sc-24988, as diluent), Western Blotting Luminol Reagent: sc-2048 and Lambda Phosphatase: sc-200312A. 2) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.